Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

ENGLISH

CAUTION - PLEASE READ CAREFULLY

Argen products should only be used by dental specialists with experience in maxillary implantology and other specialities, such as dental diagnosis, planning, dental surgery or prosthetic techniques. If in doubt regarding the product's use, please contact the manufacturer. Argen products designed for single-use only must never be re-used. If reused, there is a risk that product damage and deterioration of its characteristics could lead to prosthetic solution failure and / or other deterioration of the patient's health such as patient tissue infection.

All Argen components must be dry fitted before use to check that the fit is correct. The clinician is responsible for the correct application of restorative products as planning and procedures are under his/her control. This is the reason why only dental specialists with the appropriate experience and training should work with these products. In case of doubt please contact the manufacturer or your local distributor.

U.S. Federal law restricts those devices to sale by or on the order of a licensed dentist / physician.

We recommend an annual inspection of the prosthetic restoration by the dentist and the laboratory. This yearly inspection must consist of a screw check. If the screws are subject to unusual wear, the complete integrity of the implant abutment should be checked. New screws must be used for any adjustment correction, or replacement. Failure to follow this instruction puts the patient at risk and will void the warranty.

During any intraoral use and manipulation all Argen products must be secured to prevent aspiration due to its shape and size.

INDICATIONS

Implant abutments are used for prosthetic restorations of dental implants or for assisting procedures in the dental laboratory.

CONTRAINDICATIONS

All materials used are biocompatible; however, some patients may be allergic or hypersensitive to any of the materials and their components (specified in the table).

COMPATIBILITY INFORMATION

All Argen components are available for a variety of connections. For compatibility with dental implants and implant analogues please follow our catalogue or guides or contact your local distributor.

STERILISATION AND RE-USE

All products are supplied NON STERILE. To sterilise, we recommend autoclaving the product at 121°C for 30 minutes and letting dry for 30 minutes (in accordance with the UNE-EN ISO 17665-1:2007 standard). Some devices are marked for "Single use only" because it is difficult or impossible to clean and decontaminate a used device, reuse can lead to cross-infection. Furthermore, any attempt to reuse a device greatly increases the risk of mechanical failure caused by material fatigue. Any warranty claim resulting from the reuse of a single-use device will not be accepted.

SHELF LIFE

- ABUTMENTS/TIBASES:

The abutments and Ti-Base are intended to be handled and manipulated by professional dental technicians before being shaped into their final form. They are labelled "for single use", but do not degrade over time and are not susceptible to variations in environmental or storage conditions. Therefore, neither the shelf life of the device nor the shelf life of the package is limited and no "use by" date is provided on the label.

- SCREWS:

Clinical screws are intended to be handled by professionals for implant-prosthesis fixation, and once assembled, they are not susceptible to variations in environmental conditions. Their shelf life is not limited. Therefore, a date is not labelled.

ANALOGUE / Class I	Surgical stainless steel AISI-303	In the autoclave	Not recommended	⊗
TRANSFER / Class I	Body: surgical stainless steel AISI-303 screw: Titanium alloy ELI Ti-6Al-4V	In the autoclave	After sterilisation	⊗
SCREW / Class IIb	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
Ti BASE / Class IIb	Titanium alloy ELI Ti-6Al-4V	Autoclave before use on the patient	Single use	⊗
TEMPORARY ABUTMENT*/Class IIb	Titanium alloy ELI Ti-6Al-4V Polyether ether ketone (PEEK)	Autoclave before use on the patient	Single use	⊗
TOOLS / Class I	Handles: surgical stainless steel SUS 316L Tips: stainless steel AISI 420 MOD	In the autoclave	In case of intraoral use: after sterilisation	
TOOLS/Class IIa	Screwdrivers: surgical stainless steel SUS 316L	In the autoclave	In case of intraoral use: after sterilisation	

The Argen Corporation

8515 Miralani Drive, San Diego, CA 92126
argen.com | customercare@argen.com | (800) 255-5524

L03899 · Ver. 1 01/22

Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

ENGLISH

APPLICATIONS AND DIRECTIONS FOR USE – ANALOGUES

The analogues are for simulating the dental implant connection and position on a stone master model. Before any laboratory manipulation ensure that the anti-rotation and anti-removal parts on the analogue are securely affixed. Check that the analogue and the prosthetic element connections match in size and type before tightening. A sealed and passive connection is needed. Do not re-use an analogue, because such a connection may be impossible to achieve with a product that may no longer conform to specifications.

APPLICATIONS AND DIRECTIONS FOR USE – TRANSFER

The transfer is used in the closed and open tray technique, for the transfer of the implant position from the intraoral situation to the model in the dental laboratory. Before being used, ensure that the connection seat for the implant is clean. Any traces of dirt could affect subsequent alignment of the prosthesis. Check compatibility with the implant model to which it will be connected. Remove the abutment, clean the connection with water/air and dry. Prepare the transfer in line with the implant connection and platform size. Place it on the implant and check that it is in the correct position. Put the screw in and tighten manually.

- For the closed tray technique choose the short screw. After tightening lock the hex of the screw with wax. Use hydrocolloid, polyethylene or soft silicone. After hardening place the transfer in its groove and check the stability before sending to the laboratory. This technique is strongly recommended only for implants without axial divergence.
- For the open tray technique choose the long screw. After tightening position the impression tray such as to ensure full access to transfer screws from the outside. Before taking impressions the transfers can be ferulised with acrylic resin placed on dental floss between consecutive implants. After hardening release all screws and remove the impression tray.

APPLICATIONS AND DIRECTIONS FOR USE – CLINICAL SCREW

The screw is for fixing prosthetic restorations and auxiliary abutments over the implant or analogue. Make sure to secure the parts with their corresponding screws and observe the specified torque value placed on the label or ask your local Argen distributor. For best results, the following conditions must be adhered to meticulously:

- Use the suitable model key and size for tightening and unscrewing. If in doubt, check that the next size key does not fit into the seat. The screwdriver should be placed on the longitudinal axis of the prosthesis/implant assembly. A new screw should be used when assembling a prosthesis for the first time and for every check thereafter.
- For immediate load prostheses screw manually, avoiding excessive torque, and prevent the implant from turning while screwing.
- When transferring to the patient, do not use the same screw that was used in the laboratory.
- Make sure the correct model of screw is used for each case.
- **Positioning the patient to avoid aspiration should be avoided during screwing/unscrewing.**
- Check compatibility of the screw with the implant model to which it will be connected.

Some screw models are available with DLC coating [Diamond-Like Carbon, Chromium Carbonitride coating – CrCN] to provide a low friction surface, which improves the preload of the screw and provides better hold. The Chromium Carbonitride coating is non-toxic and has excellent mechanical, tribological and biological properties. DLC is often used in implants that make contact with blood such as stents and heart valves and to reduce wear in load bearing joints.

APPLICATIONS AND DIRECTIONS FOR USE – TI-BASE

The TiBase is used for prosthetic restorations prepared by dental technicians in a dental laboratory. The main use of the titanium TiBase is to support the bridge or individual zirconium dioxide elements made by CAD/CAM systems or with manual milling machines. It can also be used as a connector between full-shape screwed crowns (from any material) and an implant, as a narrow abutment in some cases. **Contraindications:** do not use the titanium TiBase for restorations with cantilever on a single implant, bruxism, insufficient space or direct metal-to-TiBase casting.

To improve cement adhesion we recommend that it be thoroughly cleaned and degreased before cementation. You can use any implant luting cement (duo or chemical cured) extra-orally following the recommendations of the cement manufacturer. The internal surface of the zirconium structure (cementation zone) should be sandblasted and cleaned/degreased. For secure grip the diameter and height of the TiBase should not be reduced.

APPLICATIONS AND DIRECTIONS FOR USE – MULTI-UNIT ABUTMENT

These abutments are used for prosthetic restorations prepared by dentist in dental clinic. They are a pre-manufactured dental implant abutment directly connected to the dental implant intended for use as an aid in prosthetic rehabilitation for multiple fixed-removable or screw-retained prostheses. We strongly recommend checking the final abutment seating using radiographic imaging.

In the case of using the straight MULTI-UNIT ABUTMENT: place the appropriate abutment. Use plastic holder to facilitate the insertion. Tighten the abutment using manual torque wrench (for recommended torque use the Clinical Screw Recommended Torque document). In case of using the angled MULTI-UNIT ABUTMENT: place the appropriate abutment using metal holder at desired position and tighten the abutment using the system specific screwdriver in the manual torque wrench (for recommended torque use the Clinical Screw Recommended Torque document). Take an impression using the standard procedure and provisionalise the patient. If temporary prosthesis is not necessary, place the healing caps. In using the UniAbutment®: place the adequate abutment choosing the correct type (20° or 45° - depending on implant divergence) and height. Use the attached holder screw to facilitate insertion. Take impression using standard techniques (for recommended torque use the Clinical Screw Recommended Torque document).

The Argen Corporation

8515 Miralani Drive, San Diego, CA 92126
argen.com | customercare@argen.com | (800) 255-5524

L03899 · Ver. 1 01/22

Prosthetic Fixtures for Dental Implants

INSTRUCTIONS FOR USE

ENGLISH

APPLICATIONS AND DIRECTIONS FOR USE – TEMPORARY ABUTMENT

The temporary abutment is used for prosthetic restorations prepared by dentists at the chairside, it works as support for provisional screwed crowns, bridges and prosthesis. Make sure that the engaging parts of the abutment are correctly aligned with the implant retention parts. We strongly recommend using radiography to verify the correct settlement after tightening. **Warning:** Any grinding should be performed outside the mouth. If this is not possible, adequate suction must be used. Positioned the patient so that debris is not inhaled or swallowed. Use carbide burs and carborundum discs with water cooling. Before placing the temporary restoration check the torque according to implant connection and size. It is important to seal the access hole with Teflon tape or gutta-percha and composite. Any material excess near the implant platform should be removed, as it can provoke peri-implantitis and can lead to implant loss. It is very important to check static and dynamic occlusion avoiding excessive loading forces which may prolong or prevent the proper osseointegration, especially in cases of immediate load.

APPLICATIONS AND DIRECTIONS FOR USE – ANGLED BASE

The angled Ti-Base is used for prosthetic restorations prepared with axially tilted implants. The main use of the titanium Ti-Base is to correct this axis to support the bridge or individual zirconium dioxide elements fabricated in CAD/CAM or on the manual milling-copying machine. It can also be used as a connector between screw-retained full crowns (of any material) and an implant, or in some cases as a narrow clinical abutment. Contraindications: Do not use the titanium interface for single-unit restorations with extensions on a single implant, patients with bruxism, in cases of insufficient vertical space or for metal overcasts.

To improve the adhesion of the cement it is recommended to clean and degrease the surfaces prior to cementation. Any implant cement (dual-curing or self-adhesive) can be used for cementing on the model according to the cement manufacturer's instructions. The internal surface of the zirconium structure (cementation zone) area should be sandblasted and cleaned/degreased. For secure bonding, reduction of the interface diameter and height is not recommended.

EXPLANATION OF THE SYMBOLS THAT APPEAR ON LABELS

 For single use only	LOT Batch No.	REF Reference
 NON Sterile product	 Read instructions	CE0051
Rx Only Use only under prescription	Medical product compliance marking controlled by the notified body number 0051 (Class IIb)	
 Manufacturer information	Manufacturer: Terrats Medical SL Calle Mogoda 75-99 08210 – Barberá del Vallés - Barcelona (Spain) Tel + 34 935 646 006. Fax +34 935 647 317 info@dess-abutments.com	

All DESS products are manufactured in accordance with ISO 9001 and ISO 13485 and Class IIb items bear the CE marking. 3i® and 3i® Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA. ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XIVE® and FRIADENT® and UniAbutment® are registered trademarks of DENTSPLY Implants BioHorizons® is a registered trademark of BioHorizons, INC. Brånemark System®, NobelReplace® and Replace Select™ and NobelActive™ and NobelReplace® CC and NobelSpeedy™ and MULTI-UNIT are registered trademarks of Nobel Biocare AB, Sweden. ITI®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland. MIS® is a registered trademark of MIS Implant Technologies, Ltd. Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA. Dyna® is a registered trademark of Dyna Dental Engineering BV. Locator® is a registered trademark of Zest Anchors INC.

The photographs of the products are for information purposes only. The type, reference and connection of each element should be checked. Some Argen items may not be available in your country. These instructions replace all previous editions. For any detailed information about Argen products please contact your local distributor.

Medical product compliance marking controlled by the notified body number 0051 (Class IIb)

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