INSTRUCTIONS FOR USE

SUPER TRANSLUCENCY ZIRCONIA

5855 Oberlin Drive
San Diego, CA 92121-4718 USA
EC REP MDSS, Schiffgraben 41
D-30175 Hannover, Germany

SINTERING ZIRCONIA

Stage	Program	Rate/Minute	Temperature	Time (Minutes)
1	Heating Ramp	7°C/Minute	900°C	128 Minutes
2	Heating Ramp	10°C/Minute	1450°C	55 Minutes
3	Heat Soak		1450°C	120 Minutes
4	Cooling Ramp	6°C/Minute	1000°C	75 Minutes
5	Cooling Ramp		Natural Cool	180 Minutes

*After this controlled cooling segment, the framework can cool naturally

Please refer to the complete SDS sheet provided with your order.

Technical Support
For further questions or technical support, please contact Technical Support staff at info@alloysonline.com

ADJUSTING

- Sandblast the interior and exterior of the restoration with aluminum oxide (50 micron grit size) at 40psi.
- Only use burs specifically designed for adjusting Zirconia. Always ensure that Zirconia is wet during the grinding process. A high-speed wet hand piece, at low speed, is recommended during the adjusting process in order to keep heat to a minimum.
- Avoid grinding the basal grooves and tooth connections.
- If possible, smooth rough or sharp edges.

MATERIAL PROPERTIES

STRENGTH

Flexural Bending Strength – Zirconia Density
≥ 6.00 g/cm³

COMPOSITION

ZrO ₂ +HfO ₂ +Y ₂ O ₃	>99 wt%	THERMAL EXPANSION COEFFICIENT	25-500°C = 10µm/m·°C
Y ₂ O ₃	8.5 - 10.0 wt%		
HfO ₂	< 5 wt%		
Al ₂ O ₃	< 0.1 wt%	Type/class	Type II/Class 4
Fe ₂ O ₃	< 0.1 wt%		BS EN ISO 6872:2015

HAZARDS IDENTIFICATION - EMERGENCY OVERVIEW

Specific Physical Form: Solid Block or Slab	Eye Contact: Mechanical eye irritation: Signs/symptoms may include pain, redness, tearing and corneal abrasion.	Inhalation: During grinding, scraping or sanding, inhalation of particles may occur, resulting in upper respiratory tract irritation. Signs/symptoms may include cough, sneezing, nasal discharge, headache, hoarseness, and nose and throat pain.
Odor, Color, Grade: White odoreless block	Skin Contact: Mechanical Skin irritation: Signs/symptoms may include abrasion, redness, pain, and itching.	Ingestion: No health effects are expected.
General Physical Form: Solid		
Immediate health, physical, and environmental hazards: No immediate health, physical, or environmental hazards are anticipated.		

INSTRUCTIONS FOR USE

(Super Translucent) dental zirconia is indicated for the production of full contour restorations. The following instructions provide general guidelines for handling, designing, milling, coloring, sintering and adjusting of material and should be followed very carefully to avoid any loss of aesthetics, fit, durability or overall quality.

INDICATIONS FOR USE

ArgenZ Anterior (super translucent) zirconia can be used for the production of full contour single unit anterior and posterior restorations, and 3-unit anterior bridges.

HANDLING

Inspect each shipment for damage and do not use damaged discs for the production of dental restorations. Store in a cool, dry, temperature-stable environment (between 5°C and 50°C) in the original packaging.

DESIGNING

Noncompliance with these guidelines could result in an unfit or failed restoration.

Design Option	Design Guidance
Drill Compensation	Drill compensation must be activated for all substructures milled from a solid structure.
Cement Gap	The distance where the coping intersects the die at the margin area. Use this setting to control margin fit.
Extra cement Gap	The distance between the coping walls and the die. Use this setting to control internal fit.
Distance to Margin Line	The distance from the margin outer line to the start of the interior wall of the coping.

Smooth Distance

The distance from the margin line to the margin engagement point. Should be set at 0.20mm.

Drill Radius

The drill radius is the size of the smallest end mill used to mill the pattern.

Drill Compensation Offset

The distance from the margin line to the area affected by drill compensation. Should be a minimum of 0.5mm.

Margin Line Offset

The effective thickness of the margin line and should not be less than 0.16mm. Thinner margin lines will result in a higher failure rate.

Offset Angle #1

The offset angle should not be less than 65°.

Extension Offset

The extension offset should not be less than 0.01mm.

Wall Thickness

A nominal wall thickness of 0.5mm will ensure a consistently quality product. Reducing this value could result in fractures or holes in the framework.

Bridge Connectors

Recommended Anterior restorations: 9mm² minimum.

MILLING

Pre-sintered (or "green") zirconia material has an inherent shrinkage rate associated with each production lot. This shrink rate, usually formatted as 1.XXXX, can be found on the side of the actual disc. This number MUST be input into the milling preparation software to ensure the accuracy of the eventual restoration.

When milling, always follow these general guidelines:

- Only use sharp end mills with carbide or diamond coating.
- Do not use any restoration that has chips and/or cracks. Remove the units from the disc using a handpiece with a diamond-coated burr.
- Smooth the support areas with a medium-grit rubber polishing wheel.
- Remove any residual zirconia dust with an art brush.
- If a wet mill is used make sure all the zirconia is completely dry before sintering. Air dry for at least 15 minutes prior to sintering. Damp zirconia will crack if placed in the sintering oven.

COLORING

Zirconia is compatible with all major dental zirconia coloring systems.

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