

# ArgenZ Zirconia

## Recommended Settings Parameters

---

### Material Selection Indications

**ArgenZ ST Multilayer** – Cosmetic Zirconia  
Super translucent multilayer zirconia for the esthetic zone.

**ArgenZ Anterior** – Cosmetic Zirconia  
Super translucent zirconia for the esthetic zone.

**ArgenZ HT+ Multilayer** – Restorative Esthetic Zirconia  
High translucent plus multilayer zirconia for esthetic zone and large span bridges.

**ArgenZ HT+** – Restorative Esthetic Zirconia  
High translucent plus zirconia copings and bridge frameworks.

**ArgenZ Ultra** – Restorative Zirconia  
Frameworks where extra strength and opacity is needed.  
(i.e. implant abutments and blocking out dark shades)

### Cut-off Time

Files must be submitted by **12:00pm PST (3:00pm EST)**  
1-8 unit cases will be received within 2 business days.  
9-14 unit cases will be received within 3 business days.

### Material Selection Guide

When submitting your files, please ensure you have indicated the material type (ST Multilayer, Anterior, HT+ Multilayer, HT+, Ultra).

Zirconia Type	# of Units	Type	Strength
ArgenZ ST Multilayer	Single up to 3 units	Super Translucency	850 MPa
ArgenZ Anterior	Single up to 3 units	Super Translucency	765 MPa
ArgenZ HT+ Multilayer	Single up to 14 units	High Translucency Plus	1250 MPa
ArgenZ HT+	Single up to 14 units	High Translucency Plus	1250 MPa
ArgenZ Ultra	Single up to 14 units	High Strength	1400 MPa

**3Shape users: upload via ArgenLink.**

**All other users: go to [argen.com](https://argen.com)**

Files can be uploaded one at a time, or in batches, making the upload process easy!  
Please call Argen Digital for pricing details.

# ArgenZ Zirconia

## Recommended Settings Parameters

### Designing in 3Shape

#### Design

Wall Thickness: **0.45 mm minimum**

Margin Thickness: **0.2 mm**

Margin Angle: **65° degrees**

**NOTE:** If walls/margins are too thin, failures may occur during the milling process (i.e. cracks, chipped margins).

#### Fit Settings

Drill Comp should be checked "ON"

Drill Comp offset setting: **0.60 mm**

Drill Radius: **0.50 mm**

Cement Gap: **0.045 mm** (near margin)

Extra Cement Gap: **0.065 mm**  
(die spacer)

Distance to Margin: **1.0 mm**

**NOTE:** Cement gap settings are adjustable. Increase values to give a looser fit, decrease values to get a tighter fit.

### Designing in Exocad

#### Design

Minimal Thickness: **0.45 mm minimum**

Horizontal Crown Margin: **0.2 mm**

Angled Crown Margin: **0.05 mm**

Angle: **65° degrees**

**NOTE:** If walls/margins are too thin, failures may occur during the milling process (i.e. cracks, chipped margins).

#### Fit Settings

In the wizard under "Crown Bottoms" check the box "Anticipated Milling" and set the Diameter to: **1.2 mm**

Gap thickness of Cement: **0.045 mm**  
(near margin)

Add distance occl: **0.065 mm**  
(die spacer)

Add distance X/Y: **0.065 mm**  
(die spacer)

**NOTE:** In the Wizard Merging, make sure you check the box "Optimize for Selective Laser Melting". Cement gap settings are adjustable. Increase values to give a looser fit, decrease values to get a tighter fit.

### Designing in Dental Wings

#### Design

Minimum Thickness: **0.45 mm**

Margin Thickness: **0.2 mm**

Emergence Angle: **65° degrees**

**NOTE:** If walls/margins are too thin, failures may occur during the milling process (i.e. cracks, chipped margins).

#### Fit Settings

Apply Tool Compensation box should be checked

Elliptic Vertical: **0.6 mm**

Elliptic Radial: **0.6 mm**

Additional Vertical: **0.1 mm**

Cement Gap: **0.045 mm**  
(near margin)

Extra Vertical Gap: **0.065 mm**  
(die spacer)

Extra Horizontal Gap: **0.065 mm**  
(die spacer)