

# 3Shape Scanners

Lab Scanners for Every Need

---



## E2 Increase productivity and texture scanning

- 2x5 MP cameras
- Accuracy (ISO 12836): 10  $\mu$ m
- Scan speed (full arch): 24 sec
- Scan Speed (full arch impression): 72 sec
- Texture: B/W
- Scanning Strategy: Standard



## E3 High performance and implant bar accuracy

- 2x5 MP cameras
- Accuracy (ISO 12836): 7  $\mu$ m
- Scan speed (full arch): 18 sec
- Scan Speed (full arch impression): 64 sec
- Texture: Color
- Scanning Strategy: Standard



## E4 Ultra-high speed and accuracy

- 4 x 5 MP cameras
- Accuracy (ISO 12836): 4  $\mu$ m
- Scan speed (full arch): 9 sec
- Scan Speed (full arch impression): 45 sec
- Texture: Color
- Scanning Strategy: Die-in-model



## F8 Dual model scanning

- Dual model scanning in one workflow
- 4 x 5 MP cameras
- Accuracy (ISO 12836): 4  $\mu$ m
- Scan speed (full arch): 9 sec
- Scan Speed (full arch impression): 45 sec
- Texture: Color
- Scanning Strategy: Die-in-model