

ARGEN IS IMPLANT SYSTEMS



PRODUCT COMPATIBILITY



argen.com | (800) 255-5524

ArgenIS

Implant System and Scan Bodies Compatibility

Implant Info			Indication Compatibility				Scan Body Compatibility						
System	Code	Size	Ti	Hybrid	ASC	Multi-Unit Abutment	ArgenIS **	Multi-Unit ***	NT **	3Shape	Elos	Elos IO	Dess IO
Astra Tech EV	SEV	3.0, 5.4	X	X			X		X	X		X	X
	SEV	3.6, 4.2, 4.8	X	X	X		X		X	X		X	X
Astra Tech OsseoSpeed® TX*	AO	3.0 NP		X									X
	AO	3.5/4.0 RP, 4.5/5.0 WP	X	X	X	X	X	X	X	X	X	X	X
BioHorizons®*	BIO	3.0	X				X		X	X			X
	BIO	3.5 NP, 4.5 RP, 5.7 WP	X	X			X		X	X	X	X	X
Biomet 3i Certain®*	BC	3.4 NP, 4.1 RP, 5.0 WP	X	X	X		X		X	X	X	X	X
Biomet 3i Osseotite®*	BO	3.4 NP, 5.0 WP		X	X				X	X	X	X	X
	BO	4.1 RP	X	X	X				X	X	X	X	X
CAMLOG®*	CAM	3.3, 5.0	X	X			X		X	X	X	X	X
	CAM	3.8, 4.3	X	X	X		X		X	X	X	X	X
Dentsply Friadent XiVE®*	DX	3.4 NP	X		X				X	X	X		X
	DX	3.8 RP, 4.5 WP	X		X				X	X	X	X	X
	DX	5.5	X						X	X	X	X	X
Hiossen ET/HG*	HIO	3.5, 4.0/4.5/5.0/6.0	X	X			X		X				X
MegaGen AnyRidge®*	MA	Universal	X	X	X		X						X
Neodent® Grand Morse™*	NGM	Universal	X	X	X	X	X	X					X
NobelActive®*	NC	3.5 NP, 4.3/5.0 RP	X	X	X	X	X	X	X	X	X	X	X
	NC	5.5 WP	X	X			X			X			X
Nobel Brånemark System®*	NB	3.5 NP, 4.1 RP	X	X	X				X	X	X	X	X
	NB	5.1 WP	X	X					X	X	X	X	X
NobelReplace®*	NR	3.5 NP, 4.3 RP, 5.0 WP	X	X	X	X	X	X	X	X	X	X	X
	NR	6.0		X			X		X	X			X
Straumann® Bone Level*	SB	3.3 NC, 4.1/4.8 RC	X	X	X	X	X	X	X	X	X	X	X
Straumann® synOcta®*	SS	3.5 NN		X					X	X			X
	SS	4.8 RN, 6.5 WN	X	X	X		X		X	X	X	X	X
Zimmer Tapered Screw-Vent®*	ZS	3.5 NP, 4.5 RP	X	X	X	X	X	X	X	X	X	X	X
	ZS	5.7 WP	X	X	X		X		X	X	X	X	X

* Not trademarks of The Argen Corporation. ** Sold by Argen. *** Sold by Argen, one size fits all.

MKT-18982-Ver4-09/24

Ti Bases Gingival Heights

All platforms available in engaging and non-engaging

Series	Size	Gingival Heights
Astra Tech EV* (SEV)	3.6	1.0mm, 3.0mm
	4.2	1.0mm, 3.0mm
	4.8	1.0mm, 3.0mm
Astra Tech OsseoSpeed® TX* (AO)	3.0 NP	1.0mm, 3.0mm
	3.5/4.0 RP	1.1mm, 3.0mm
	4.5/5.0 WP	0.9mm, 3.0mm
BioHorizons®* (BIO)	3.5 NP	0.8mm
	4.5 RP	0.3mm
	5.7 WP	0.3mm
Biomet 3i Certain®* (BC)	3.4 NP	0.3mm, 1.5mm
	4.1 RP	0.3mm, 1.5mm
	5.0 WP	0.3mm, 1.5mm
Biomet 3i Osseotite®* (BO)	3.4 NP	0.3mm
	4.1 RP	0.3mm
	5.0 WP	0.3mm
CAMLOG®* (CAM)	3.3	0.5mm
	3.8	0.3mm, 1.5mm
	4.3	0.3mm, 1.5mm
	5.0	0.3mm, 1.5mm
Hiossen ET/HG* (HIO)	3.5	1.0mm, 3.0mm
	4.0-6.0	1.0mm, 3.0mm
MegaGen AnyRidge®* (MA)	Universal	2.0mm
Neodent® Grand Morse™* (NGM)	Universal	0.8mm, 1.5mm
NobelActive®* (NC)	3.5 NP	1.2mm, 3.0mm
	4.3/5.0 RP	1.2mm, 3.0mm
	5.5 WP	1.5mm

Series	Size	Gingival Heights
Nobel Brånemark System®* (NB)	3.5 NP	0.3mm
	4.1 RP	0.3mm
	5.1 WP	0.3mm
NobelReplace® Select™* (NR)	3.5 NP	0.3mm
	4.3 RP	0.3mm
	5.0 WP	0.3mm
Straumann®* Bone Level (SB)	3.3 NC	1.2mm, 3.0mm
	4.1/4.8 RC	1.2mm, 3.0mm
Straumann® synOcta®* (SS)	3.5 NN	0.5mm
	4.8 RN	0.4mm
	6.5 WN	0.3mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5 NP	0.8mm, 1.5mm
	4.5 RP	0.3mm, 1.5mm
	5.7 WP	0.3mm, 1.5mm

MKT-18953-Ver4-09/24

Multi-Unit Abutment

Gingival Heights and Angles

Series	Size	Multi-Unit Platform	Angle	Gingival Heights
Astra Tech OsseoSpeed® TX* (AO)	3.5/4.0	4.8	Straight	1.5 mm, 2.5 mm
	4.5/5.0	4.8	Straight	1.5 mm, 2.5 mm
Neodent® Grand Morse™*	Universal	4.8	Straight	1.5 mm, 2.5 mm
	Universal	4.8	17°	2.5 mm, 3.5 mm
	Universal	4.8	30°	3.5 mm, 4.5 mm
NobelActive®* (NC)	3.5 NP	4.8	Straight	1.5 mm, 2.5 mm
	3.5 NP	4.8	17°	2.5 mm, 3.5 mm
	3.5 NP	4.8	30°	3.5 mm, 4.5 mm
	4.3/5.0 RP	4.8	Straight	1.5 mm, 2.5 mm
	4.3/5.0 RP	4.8	17°	2.5 mm, 3.5 mm
	4.3/5.0 RP	4.8	30°	3.5 mm, 4.5 mm
NobelReplace® Select™* (NR)	3.5 NP	4.8	Straight	1.0 mm, 2.0 mm
	4.3 RP	4.8	Straight	1.0 mm, 2.0 mm
	5.0 WP	4.8	Straight	1.0 mm, 2.0 mm
Straumann® Bone Level* (SB)	3.3 NC	4.8	Straight	1.5 mm, 2.5 mm
	3.3 NC	4.8	17°	2.5 mm, 3.5 mm
	3.3 NC	4.8	30°	3.5 mm, 4.5 mm
	4.1/4.8 RC	4.8	Straight	1.5 mm, 2.5 mm
	4.1/4.8 RC	4.8	17°	2.5 mm, 3.5 mm
	4.1/4.8 RC	4.8	30°	3.5 mm, 4.5 mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5	4.8	Straight	1.0 mm, 2.0 mm
	3.5	4.8	17°	2.5 mm, 3.5 mm
	3.5	4.8	30°	3.5 mm, 4.5 mm
	4.5	4.8	Straight	1.0 mm, 2.0 mm
	4.5	4.8	17°	2.5 mm, 3.5 mm
	4.5	4.8	30°	3.5 mm, 4.5 mm

MKT-18965-Ver4-09/24