


ARGEN IS



IMPLANT SYSTEMS

Quality implant parts from a partner you trust.



argen.com | (800) 255-5524

Table of Contents

Why ArgenIS?	4-10
---------------------------	-------------

Implant System Components	11-33
--	--------------

Astra Tech EV*	11
Astra Tech OsseoSpeed TX®*	13
BioHorizons®*	15
Biomet 3i Certain®*	16
Biomet 3i Osseotite®*	17
CAMLOG®*	18
Dentsply Friadent XiVE®*	20
Hiossen ET/HG*	21
MegaGen AnyRidge®*	22
Neodent® Grand Morse™*	23
NobelActive®*	24
NobelReplace Select®*	26
Nobel Brånemark System®*	28
Straumann® Bone Level*	29
Straumann® synOcta® Tissue Level*	31
Zimmer Tapered Screw-Vent®*	32

Multi-Unit Abutments (MUAs)	34-40
--	--------------

Universal Compatibility	34
Astra Tech OsseoSpeed TX®*	35
Neodent® Grand Morse™*	36
NobelActive®*	37
NobelReplace® Select™*	38
Straumann® Bone Level*	39
Zimmer Tapered Screw-Vent®*	40

*Not a trademark of The Argen Corporation.



Tools 42-43

Torque Wrenches 42

Adapters..... 42

ArgenIS Screwdrivers..... 42

Other Compatible Screwdrivers..... 43

Screwdriver Handles and Inserts 43

Multi-Unit Abutment Keys..... 43

Compatibility Charts 44-46

ArgenIS Implant System and Scan Bodies 44

Ti Bases Gingival Heights 45

Multi-Unit Abutment Gingival Heights and Angles 46



ArgenIS custom abutments are manufactured with one thing in mind: **the patient.**

Argen's advanced engineering and state-of-the-art manufacturing produces precise, industry leading tolerances between the implant and abutment to ensure excellent fit every time. Argen's custom abutments are compatible with the leading implant systems for maximum efficiency and function. We provide extensive warranty offerings for increased customer satisfaction and your complete peace of mind.

Argen's digital technology services help your lab elevate its capabilities by delivering a digital workflow that is unrivaled in the industry. We can implement a completely custom workflow based on your needs, from scan to model, design of the abutment, to the coping and crown. Experience best-in-class speed with our fast, two-day turnaround time.

At Argen, we understand that dental laboratories are building a smile, not just abutments, one patient at a time.

Argen's Implant Specific Abutments

Argen's custom abutments and screws are certified OEM-equivalent; they function comparably to OEM (Original Equipment Manufacturer) parts. Our manufacturing processes and standards are designed to create the highest quality implant solutions compatible with the leading implant systems.

Laboratories, doctors, and patients can have complete peace of mind knowing the performance of Argen's implant parts is assured and covered by an extensive warranty.

Why dental laboratories choose Argen's Implant Specific Abutments

	Argen Implant Specific Abutments	Various OEM Abutments
FDA 510(k) Cleared	✓	✓
MDSAP Certified	✓	Varies
ISO 13485	✓	Varies
Warranty (Implant & Abutment)	✓	Varies
Strict Tolerances	✓	✓
True-to-Design Margins	✓	✓
Split Files	✓	Varies
FDA 510(k) Cleared 2 Final Screws	✓	Varies
Better Value	✓	---



Extensive Warranty

Custom Abutment and Bar Warranty

Argen provides a limited warranty for all Argen manufactured custom abutments and specified bars. In addition, if an implant supplier does not honor their warranty due to the use of an Argen abutment or bar, Argen will replace both the abutment and the implant – certain limitations apply. Please visit argen.com for complete terms and conditions of the warranty.



Visit argen.com for complete terms of the warranty.



FDA 510(k) Cleared / ISO 13485 MDSAP

FDA 510(k) Cleared, ISO 13485 Certified, MDSAP Certified

While most OEM manufacturers have 510(k) clearance, Argen has gone a step further to become MDSAP compliant. MDSAP is a global approach to auditing and monitoring the safe manufacturing of medical devices. It is not a one-time pass or fail; Argen is routinely monitored to ensure we meet these standards. This compliance provides our customers with peace of mind knowing Argen consistently offers quality products.



Robotic Automation

Robotic Automation

Argen's custom-built 5-axis FANUC Robodrills, 6-axis FANUC robots, and software work together to streamline our abutment milling process. The integrated solution automates our abutment manufacturing, ensuring consistency and quality throughout the process.



Advanced Milling Strategies

Advanced Milling Strategy

- Advanced milling strategies for true-to-design margins and finishes
- Automated nesting and strategy selection for optimal efficiency
- Rapid toolpath calculation and machine motion optimization for increased speed to completion
- Post-machining probing for accurate inspection



End-to-End Manufacturing

Vertical Integration

Argen's engineered end-to-end system provides complete oversight of the production of blanks, milled custom abutments, and final screws. We even manufacture the tooling used in the milling process. This level of control allows us to produce the most precise and consistent OEM-equivalent parts.

Argen – A Partner You Can Trust

Argen is committed to doing whatever it takes to make dental labs successful. Our extensive abutment warranty includes the abutment and implant.

“ Before using Argen custom abutments, we struggled to design abutments to our exact specifications due to manufacturing limitations. Argen solved this by allowing us to create precise designs with minimal hand customization.

We value their ease of uploads, exceptional customer service, accurate design-to-abutment match, high-quality finish, and the inclusion of two screws. Their fast turnaround, consistent quality, and variety of implant systems make them a trusted partner.

I would recommend Argen for their unmatched quality, consistency, and value — they've transformed our workflow and results.


”

Brad T Alioto
CAD/CAM Manager, Owner
Expertec Dental Laboratory, Inc.



Watch Argen's ISA manufacturing video.
argen.com/custom-abutments

Designed Together. Manufactured Together. Shipped Together.



Crown

Complete your case efficiently in one digital workflow using an ArgenZ™ zirconia crown. Argen manufactures a variety of top-quality zirconia options to meet the esthetic demands of each case.

Abutment

Advanced engineering, state-of-the-art manufacturing and our FDA-cleared proprietary manufacturing processes ensure precise, industry-leading tolerances for an optimal fit every time.

DIM Analog

The ArgenIS Digital Implant Model (DIM) analog is used in combination with printed models in digital workflows. Consisting of two pieces, this implant analog ensures accurate alignment of its mating parts and is reusable between models of the same system and size.


Digital Model

Argen is an industry leader in digital implant models using state-of-the-art technology in 3D printing.

ArgenIS

Product Offering

Argen gives you the flexibility to scan and design your own abutment and submit your design file to argen.com for precision milling by our professional CAD/CAM engineers. Our custom abutments are milled with a 5-axis precision mill for the highest quality units, with no finishing or polishing required.



Custom Abutment Type	Order Cut-Off Time	Turnaround Time
Hybrid, Angled Screw Channel, Multi-Unit Zirconia	12 pm PST (3 pm EST)	2 Business Days <i>Zirconia bridges need 1 additional day for sintering</i>
Hybrid, Angled Screw Channel, Multi-Unit Non-Zirconia	3 pm PST (6 pm EST)	2 Business Days
Titanium Abutments	3 pm PST (6 pm EST)	2 Business Days
Titanium, CoCr Bars	3 pm PST (6 pm EST)	2 Business Days


ArgenIS Abutment Components

ArgenIS system includes: base set, lab analogs, abutment screws, and digital implant model (DIM) analogs to support leading implant systems. All lab analogs are stainless steel. Abutment screws come with recommended screw torque settings per compatible system.



Compatible with the Leading Implant Systems

- ✓ Astra Tech EV*
- ✓ Astra Tech OsseoSpeed®*
- ✓ BioHorizons®*
- ✓ Biomet 3i Certain®*
- ✓ Biomet 3i Osseotite®*
- ✓ CAMLOG®*
- ✓ Dentsply Friadent XiVE™*
- ✓ Hiossen ET / HT®*
- ✓ MegaGen AnyRidge®*
- ✓ Neodent® Grand Morse™*
- ✓ NobelActive®*
- ✓ NobelReplace™*
- ✓ Nobel Brånemark®*
- ✓ Straumann® Bone Level*
- ✓ Straumann® synOcta® Tissue Level*
- ✓ Zimmer Tapered Screw-Vent®*








*Not a trademark of The Argen Corporation.

ArgenIS

Compatible with **ASTRA TECH EV***

SEV Series




Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
	3.0	Single	136551	132380
	3.6	Single	136552	136304
	4.2	Single	136553	136305
	4.8	Single	136554	136306
	5.4	Single	136555	136307
				1 each
				136719
				2 each
				136718



Straight Ti Bases - Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU	
	3.6	1.0mm	E	131344
	3.6	1.0mm	NE	131341
	3.6	3.0mm	E	131348
	3.6	3.0mm	NE	131347
	4.2	1.0mm	E	131345
	4.2	1.0mm	NE	131342
	4.2	3.0mm	E	131350
	4.2	3.0mm	NE	131349
	4.8	1.0mm	E	131346
	4.8	1.0mm	NE	131343
	4.8	3.0mm	E	131352
	4.8	3.0mm	NE	131351



Lab Analogs

Size	SKU	
	3.0	136508
	3.6	131338
	4.2	131339
	4.8	131340
	5.4	136509



DIM Analogs

Size	SKU	
	3.0	121018
	3.6	121019
	4.2	121020
	4.8	121021
	5.4	121022

Continued on next page

*Not a trademark of The Argen Corporation.







ArgenIS

Compatible with **ASTRA TECH EV***

SEV Series






Angled Screw Channel Ti Bases – Comes with screw

	Size	Engaging/Non-Engaging	SKU
	3.6	E	127500
	3.6	NE	127501
	4.2	E	127502
	4.2	NE	127503
	4.8	E	127504
	4.8	NE	127505






Angled Screw Channel Replacement Screw

	Size	SKU
	3.6	127506
	4.2	127507
	4.8	127508



Ti Base and Ti Custom Abutment Replacement Screw

	Size	SKU
	3.0	132693
	3.6	131335
	4.2	131336
	4.8 / 5.4	131337

Tools



Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with ASTRA TECH OsseoSpeed TX^{®*}









AO Series




Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
 3.5/4.0 RP	Single	136526	120821	<u>1 each</u> 120866
 4.5/5.0 WP	Single	136527	120822	<u>2 each</u> 120855






Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
 3.0 NP	1.0 mm	E	131368
 3.0 NP	1.0 mm	NE	131365
 3.5/4.0 RP	1.1 mm	E	131369
 3.5/4.0 RP	1.1 mm	NE	131366
 3.5/4.0 RP	3.0 mm	E	131373
 3.5/4.0 RP	3.0 mm	NE	131371
 4.5/5.0 WP	0.9 mm	E	131370
 4.5/5.0 WP	0.9 mm	NE	131367
 4.5/5.0 WP	3.0 mm	E	131374
 4.5/5.0 WP	3.0 mm	NE	131372






Lab Analogs

Size	SKU
 3.0 NP	131362
 3.5/4.0 RP	131363
 4.5/5.0 WP	131364







DIM Analogs

Size	SKU
 3.0 NP	136276
 3.5/4.0 RP	115943
 4.5/5.0 WP	115944





Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/ Non-Engaging	SKU
 3.5/4.0 RP	E	127509
 3.5/4.0 RP	NE	127510
 4.5/5.0 WP	E	127511
 4.5/5.0 WP	NE	127512





Angled Screw Channel Replacement Screw

Size	SKU
 3.5/4.0 RP	127506
 4.5/5.0 WP	127513



Transfer Posts – Comes with two screws

Size	SKU
 3.5/4.0 RP	131375
 4.5/5.0 WP	131376

Impression transfer packaged with open and closed tray screw.





Continued on next page

*Not a trademark of The Argon Corporation.

ArgenIS compatible with ASTRA TECH OsseoSpeed TX^{®*} AO Series






Temporary Abutments *-Comes with screw*

	Size	Engaging/Non-Engaging	SKU
	3.5/4.0 RP	E	131379
	3.5/4.0 RP	NE	131377
	4.5/5.0 WP	E	131380
	4.5/5.0 WP	NE	131378





Ti Base Replacement Screw

	Size	SKU
	3.0 NP	131359
	3.5/4.0 RP	131360
	4.5/5.0 WP	131361



Ti Custom Abutment Replacement Screw

	Size	SKU
	3.5/4.0 RP	120078
	4.5/5.0 WP	120079

Tools





Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgentIS compatible with BioHorizons®*







BIO Series




Scan Bodies		Ti	PEEK	PEEK Kit	
Size	Quantity	SKU	SKU	SKU	
	3.0	Single	136530	132381	1 each 136721
	3.5 NP	Single	136531	136308	
	4.5 RP	Single	136532	136309	2 each 136720
	5.7 WP	Single	136533	136310	







Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU	
	3.5 NP	0.8 mm	E	131634
	3.5 NP	0.8 mm	NE	131631
	4.5 RP	0.3 mm	E	131635
	4.5 RP	0.3 mm	NE	131632
	5.7 WP	0.3 mm	E	131636
	5.7 WP	0.3 mm	NE	131633







Lab Analogs

Size	SKU	
	3.0	135704
	3.5 NP	131628
	4.5 RP	131629
	5.7 WP	131630



DIM Analogs

Size	SKU	
	3.0	136318
	3.5 NP	136315
	4.5 RP	136316
	5.7 WP	136317



Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
Universal	131627

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argon Corporation.

ArgenIS compatible with Biomet 3i Certain®*

BC Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
3.4 NP	Single	136528	120823	1 each 120860
4.1 RP/5.0 WP	Single	136529	120824	2 each 120849




Ti Bases – Comes with screw

Size	Gingival Height	Engaging/Non-Engaging	SKU
3.4 NP	0.3 mm	E	131315
3.4 NP	0.3 mm	NE	131312
3.4 NP	1.5 mm	E	131319
3.4 NP	1.5 mm	NE	131318
4.1 RP	0.3 mm	E	131316
4.1 RP	0.3 mm	NE	131313
4.1 RP	1.5 mm	E	131321
4.1 RP	1.5 mm	NE	131322
5.0 WP	0.3 mm	E	131317
5.0 WP	0.3 mm	NE	131314
5.0 WP	1.5 mm	E	131323
5.0 WP	1.5 mm	NE	131320



Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
Universal	131306




Lab Analogs

Size	SKU
3.4 NP	131307
4.1 RP	131308
5.0 WP	131309



DIM Analogs

Size	SKU
3.4 NP	136265
4.1 RP	136266
5.0 WP	136267



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.4 NP	E	127486
3.4 NP	NE	127487
4.1 RP	E	127488
4.1 RP	NE	127489
5.0 WP	E	127490
5.0 WP	NE	127491



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127492

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgentIS compatible with Biomet 3i Osseotite®*

BO Series



Scan Bodies

Size	Quantity	Ti SKU
3.4 NP	Single	135140
4.1 RP	Single	135141
5.0 WP	Single	135142



Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
3.4 NP	0.3 mm	E	131332
4.1 RP	0.3 mm	E	131333
5.0 WP	0.3 mm	E	131334



Ti Base and Custom Abutment Replacement Screw

Size	SKU
Universal	131324



Lab Analogs

Size	SKU
3.4 NP	131745
4.1 RP	131325
5.0 WP	131326



DIM Analogs

Size	SKU
3.4 NP	136268
4.1 RP	136269
5.0 WP	136270



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.4 NP	E	127493
3.4 NP	NE	127494
4.1 RP	E	127495
4.1 RP	NE	127496
5.0 WP	E	127497
5.0 WP	NE	127498



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127499

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with

CAMLOG®*

CAM Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
3.3	Single	132382	136534	1 each
3.8	Single	136311	136535	136723
4.3	Single	136312	136536	2 each
5.0	Single	136313	136537	136722



Ti Bases – Comes with screw

Size	Gingival Height	Engaging/Non-Engaging	SKU
3.3	0.5 mm	E	131395
3.3	0.5 mm	NE	131391
3.8	0.3 mm	E	131396
3.8	0.3 mm	NE	131392
3.8	1.5 mm	E	131401
3.8	1.5 mm	NE	131399
4.3	0.3 mm	E	131397
4.3	0.3 mm	NE	131393
4.3	1.5 mm	E	131402
4.3	1.5 mm	NE	131400
5.0	0.3 mm	E	131398
5.0	0.3 mm	NE	131394



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.8	E	127514
3.8	NE	127515
4.3	E	127516
4.3	NE	127517



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127518



Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
3.3/3.8/4.3	131385
5.0	131386



Transfer Posts – Comes with two screws

Size	SKU
3.8	131403
4.3	131404

Impression transfer packaged with open and closed tray screw.


Continued on next page

*Not a trademark of The Argon Corporation.





Argenis compatible with

CAMLOG®*

CAM Series



Lab Analogs

	Size	SKU
	3.3	131387
	3.8	131388
	4.3	131389
	5.0	131390







DIM Analogs

	Size	SKU
	3.3	136279
	3.8	136280
	4.3	136281
	5.0	136282



Temporary Abutments - Comes with screw

	Size	Engaging/Non-Engaging	SKU
	3.8	E	131405
	3.8	NE	131406
	4.3	E	131407
	4.3	NE	131408

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with Dentsply Friadent XiVE®*

DX Series



Scan Bodies

	Size	Quantity	SKU
	3.4 NP	Single	135158
	3.8 RP	Single	135159
	4.5	Single	135160
	5.5	Single	135181



Ti Custom Abutment Replacement Screw

Size	SKU
Universal	131409

Tools

Refer to Tools section of the catalog.



Lab Analogs

	Size	SKU
	3.4	131410
	3.8	131411
	4.5	131412
	5.5	131413



DIM Analogs

	Size	SKU
	3.4	136297
	3.8	136298
	4.5	136299
	5.5	136300



Angled Screw Channel Ti Bases - Comes with screw

	Size	Engaging/Non-Engaging	SKU
	3.4	E	127519
	3.4	NE	127520
	3.8	E	127521
	3.8	NE	127522
	4.5	E	127523
	4.5	NE	127524



Angled Screw Channel Replacement Screw



Size	SKU
Universal	127492

*Not a trademark of The Argen Corporation.

ArgentIS compatible with Hiossen ET/HG*









HIO Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
 3.5	Single	136538	132383	<u>1 each</u> 136727
 4.0 - 6.0	Single	136539	136314	<u>2 each</u> 136726





Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
 3.5	1.0mm	E	131666
 3.5	1.0mm	NE	131664
 3.5	3.0mm	E	131670
 3.5	3.0mm	NE	131668
 4.0 - 6.0	1.0mm	E	131667
 4.0 - 6.0	1.0mm	NE	131665
 4.0 - 6.0	3.0mm	E	131671
 4.0 - 6.0	3.0mm	NE	131669





Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
 3.5	131660
 4.0 - 6.0	131661





Lab Analogs

Size	SKU
 3.5	131662
 4.0 - 6.0	131663



DIM Analogs

Size	SKU
 3.5	136319
 4.0 - 6.0	136320

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with MegaGen AnyRidge®*

MA Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
Universal	Single	136540	132365	1 each 136729
				2 each 136728



Ti Bases - Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
Universal	2.0mm	E	131708
Universal	2.0mm	NE	131707



Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
Universal	131705




Lab Analogs

Size	SKU
Universal	131706



DIM Analogs

Size	SKU
Universal	131774



Angled Screw Channel Ti Bases - Comes with screw

Size	Engaging/Non-Engaging	SKU
Universal	E	131711
Universal	NE	131712



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127507

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with Neodent® Grand Morse™* NGM Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
Universal	Single	136544	135834	1 each 136725
				2 each 136724



Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
Universal	0.8 mm	E	133397
Universal	0.8 mm	NE	133395
Universal	1.5 mm	E	133398
Universal	1.5 mm	NE	133396



Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
Universal	133392



Lab Analogs

Size	SKU
Universal	133393



DIM Analogs

Size	SKU
Universal	133394



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
2.5 mm	E	133400
2.5 mm	NE	133399



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127518



Transfer Posts – Comes with two screws

Size	SKU
Universal	133401

Impression transfer packaged with open and closed tray screw.

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with NobelActive®*

NC Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
3.5 NP	Single	136541	120834	1 each 120859
4.3/5.0 RP	Single	136542	120835	2 each
5.5 WP	Single	136543	121488	120848




Ti Bases – Comes with screw

Size	Gingival Height	Engaging/Non-Engaging	SKU
3.5 NP	1.2 mm	E	131434
3.5 NP	3.0 mm	E	131439
3.5 NP	1.2 mm	NE	131430
3.5 NP	3.0 mm	NE	131437
4.3/5.0 RP	1.2 mm	E	131435
4.3/5.0 RP	3.0 mm	E	131440
4.3/5.0 RP	1.2 mm	NE	131431
4.3/5.0 RP	3.0 mm	NE	131438
5.5 WP	1.5 mm	E	131436
5.5 WP	1.5 mm	NE	131432



Ti Base Replacement Screw

Size	SKU
3.5 NP	131424
4.3/5.0 RP, 5.5 WP	131425



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.5 NP	E	127525
3.5 NP	NE	127526
4.3/5.0 RP	E	127527
4.3/5.0 RP	NE	127528



Angled Screw Channel Replacement Screw

Size	SKU
3.5 NP	127506
4.3/5.0 RP	127529



Ti Custom Abutment Replacement Screw

Size	SKU
3.5 NP	120074
4.3/5.0 RP, 5.5 WP	120075



Transfer Posts – Comes with two screws

Size	SKU
3.5 NP	131441
4.3/5.0 RP	131442


Impression transfer packaged with open and closed tray screw.

Continued on next page




*Not a trademark of The Argen Corporation.

ArgenIS compatible with NobelActive®*

NC Series






Lab Analogs

Size	SKU
 3.5 NP	131426
 4.3/5.0 RP	131427
 5.5 WP	131428







DIM Analogs

Size	SKU
 3.5 NP	136283
 4.3/5.0 RP	136284
 5.5 WP	131773



Temporary Abutments *-Comes with screw*

Size	Engaging/ Non-Engaging	SKU
 3.5 NP	E	131445
 3.5 NP	NE	131443
 4.3/5.0 RP	E	131446
 4.3/5.0 RP	NE	131444

Tools




Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

ArgenIS compatible with NobelReplace® Select™*

NR Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
 3.5 NP	Single	136545	120836	1 each 120858
 4.3 RP	Single	136546	120837	2 each 120847
 5.0 WP	Single	136547	120838	





Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
 3.5 NP	0.3 mm	E	131486
 3.5 NP	0.3 mm	NE	131482
 4.3 RP	0.3 mm	E	131487
 4.3 RP	0.3 mm	NE	131483
 5.0 WP	0.3 mm	E	131488
 5.0 WP	0.3 mm	NE	131484
 6.0	0.3 mm	E	131489
 6.0	0.3 mm	NE	131485









Ti Base and Ti Custom Abutment Replacement Screw

Size	SKU
 3.5 NP	131476
 4.3 RP, 5.0 WP, 6.0	131477





Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
 3.5 NP	E	127535
 3.5 NP	NE	127536
 4.3 RP	E	127537
 4.3 RP	NE	127538
 5.0 WP	E	127539
 5.0 WP	NE	127540






Angled Screw Channel Ti Base Replacement Screw

Size	SKU
 3.5 NP	127541
 4.3 RP, 5.0 WP	127542



Transfer Posts – Comes with two screws

Size	SKU
 3.5 NP	131490
 4.3 RP	131491
 5.0 WP	131492


Impression transfer packaged with open and closed tray screw.

Continued on next page




*Not a trademark of The Argen Corporation.

ArgenIS compatible with NobelReplace® Select™*

NR Series



Lab Analogs

	Size	SKU
	3.5 NP	131478
	4.3 RP	131479
	5.0 WP	131480









DIM Analogs

	Size	SKU
	3.5 NP	136288
	4.3 RP	136289
	5.0 WP	136290
	6.0	136291



Temporary Abutments - Comes with screw

	Size	Engaging/Non-Engaging	SKU
	3.5 NP	E	131493
	3.5 NP	NE	131496
	4.3 RP	E	131494
	4.3 RP	NE	131497
	5.0 WP	E	131495
	5.0 WP	NE	131498

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argon Corporation.

ArgenIS compatible with Nobel Brånemark System®*

NB Series



Scan Bodies

Size	Quantity	SKU
3.5 NP	Single	135131
4.1 RP	Single	135132
5.1 WP	Single	135133



Ti Bases

Size	Gingival Height	Engaging/ Non-Engaging	SKU
3.5 NP	0.3 mm	E	131469
4.1 RP	0.3 mm	E	131470
5.1 WP	0.3 mm	E	131471



Ti Base Replacement Screw

Size	SKU
3.5 NP	131424
4.1 RP	131461
5.1 WP	131462



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.5 NP	E	127530
3.5 NP	NE	127531
4.1 RP	E	127532
4.1 RP	NE	127533



Angled Screw Channel Replacement Screw

Size	SKU
3.5 NP	127534
4.1 RP	127499



Ti Custom Abutment Replacement Screw

Size	SKU
3.5 NP	131424
4.1 RP	131461
5.1 WP	131462



Lab Analogs

Size	SKU
3.5 NP	131463
4.1 RP	131464
5.1 WP	131465



DIM Analogs

Size	SKU
3.5 NP	136285
4.1 RP	136286
5.1 WP	136287



Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argon Corporation.







ArgenIS compatible with Straumann® Bone Level^{*} SB Series



Scan Bodies		Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU
 3.3 NC	Single	136549	120840	<u>1 each</u> 120863
 4.1/4.8 RC	Single	136550	120841	<u>2 each</u> 120852





Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
 3.3 NC	1.2 mm	E	131513
 3.3 NC	1.2 mm	NE	131511
 3.3 NC	3.0 mm	E	131517
 3.3 NC	3.0 mm	NE	131515
 4.1/4.8 RC	1.2 mm	E	131514
 4.1/4.8 RC	1.2 mm	NE	131512
 4.1/4.8 RC	3.0 mm	E	131518
 4.1/4.8 RC	3.0 mm	NE	131516





Ti Base Replacement Screw

Size	SKU
 3.3 NC	131507
 4.1/4.8 RC	131508





Lab Analogs

Size	SKU
 3.3 NC	131509
 4.1/4.8 RC	131510







DIM Analogs

Size	SKU
 3.3 NC	136292
 4.1/4.8 RC	136293



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
 3.3 NC	E	127543
 3.3 NC	NE	127544
 4.1/4.8 RC	E	127545
 4.1/4.8 RC	NE	127546



Angled Screw Channel Replacement Screw

Size	SKU
Universal	127547



Continued on next page


*Not a trademark of The Argen Corporation.

ArgenIS compatible with Straumann® Bone Level^{*} SB Series





Ti Custom Abutment Replacement Screw

Size	SKU
 3.3 NC	120076
 4.1/4.8 RC	120077







Transfer Posts - Comes with screw

Size	SKU
 3.3 NC	131527
 4.1/4.8 RC	131528

Impression transfer packaged with open and closed tray screw.



Temporary Abutments - Comes with screw

Size	Engaging/Non-Engaging	SKU
 3.3 NC	E	131525
 3.3 NC	NE	131523
 4.1/4.8 RC	E	131526
 4.1/4.8 RC	NE	131524

Tools

Refer to Tools section of the catalog.

^{*}Not a trademark of The Argen Corporation.

ArgenIS compatible with Straumann® synOcta® Tissue Level* SS Series



Scan Bodies

			Ti	PEEK	PEEK Kit
	Size	Quantity	SKU	SKU	SKU
	3.5 NN (Hybrid)	Single	136556	NA	1 each 120864
	4.8 RN	Single	136557	120842	2 each 120853
	6.5 WN	Single	136558	120843	



Ti Bases – Comes with screw

	Size	Gingival Height	Engaging/ Non-Engaging	SKU
	3.5 NN (Hybrid)	0.5 mm	E	131553
	3.5 NN (Hybrid)	0.5 mm	NE	131550
	4.8 RN	0.4 mm	E	131554
	4.8 RN	0.4 mm	NE	131551
	6.5 WN	0.3 mm	E	131555
	6.5 WN	0.3 mm	NE	131552



Ti Base and Ti Custom Abutment Replacement Screw

	Size	SKU
	3.5 NN (Hybrid)	131542
	4.8 RN/ 6.5 WN	131543




Lab Analogs

	Size	SKU
	3.5 NN (Hybrid)	131545
	4.8 RN	131546
	6.5 WN	131547



DIM Analogs

	Size	SKU
	3.5 NN (Hybrid)	136296
	4.8 RN	136294
	6.5 WN	136295



Angled Screw Channel Ti Bases – Comes with screw

	Size	Engaging/Non-Engaging	SKU
	4.8 RN	E	127548
	4.8 RN	NE	127549
	6.5 WN	E	127550
	6.5 WN	NE	127551



Angled Screw Channel Replacement Screw

	Size	SKU
	Universal	127552

Tools

Refer to Tools section of the catalog.


*Not a trademark of The Argen Corporation.

ArgenIS compatible with Zimmer Tapered Screw-Vent®*

ZS Series



Scan Bodies			Ti	PEEK	PEEK Kit
Size	Quantity	SKU	SKU	SKU	
3.5 NP	Single	136559	120844	1 each 120865	
4.5 RP	Single	136560	120845	2 each 120854	
5.7 WP	Single	136561	120846		



Ti Bases – Comes with screw

Size	Gingival Height	Engaging/ Non-Engaging	SKU
3.5 NP	0.8mm	E	131575
3.5 NP	0.8mm	NE	131573
3.5 NP	1.5mm	E	131576
3.5 NP	1.5mm	NE	131574
4.5 RP	0.3mm	E	131579
4.5 RP	0.3mm	NE	131577
4.5 RP	1.5mm	E	131580
4.5 RP	1.5mm	NE	131578
5.7 WP	0.3mm	E	131583
5.7 WP	0.3mm	NE	131581
5.7 WP	1.5mm	E	131584
5.7 WP	1.5mm	NE	131582



Ti Base Replacement Screw

Size	SKU
Universal	131569



Ti Custom Abutment Replacement Screw

Size	SKU
Universal	120069



Lab Analogs

Size	SKU
3.5 NP	131570
4.5 RP	131571
5.7 WP	131572



DIM Analogs


Size	SKU
3.5 NP	136301
4.5 RP	136302
5.7 WP	136303

Continued on next page

*Not a trademark of The Argen Corporation.

Argenis compatible with Zimmer Tapered Screw-Vent®*

ZS Series



Angled Screw Channel Ti Bases – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.5 NP	E	127553
3.5 NP	NE	127554
4.5 RP	E	127555
4.5 RP	NE	127555
5.7 WP	E	127557
5.7 WP	NE	127558



Angled Screw Channel Replacement Screw


Size	SKU
Universal	127559



Transfer Posts – Comes with screw

Size	SKU
3.5 NP	131585
4.5 RP	131586
5.7 WP	131587

Impression transfer packaged with open and closed tray screw.



Temporary Abutments – Comes with screw

Size	Engaging/Non-Engaging	SKU
3.5 NP	E	131591
3.5 NP	NE	131588
4.5 RP	E	131592
4.5 RP	NE	131589
5.7 WP	E	131593
5.7 WP	NE	131590

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

Multi-Unit Abutments

Universal Compatibility



Multi-Unit Abutment Lab Scan Body


-Comes with screw

Multi-Unit Platform	SKU
4.8	131743



Multi-Unit Abutment Lab Scan Body and Ti Base Replacement Screw

Multi-Unit Platform	SKU
4.8	131741



Multi-Unit Abutment Transfer Post

-Comes with screw

Multi-Unit Platform	SKU
4.8	131732



Multi-Unit Abutment Temporary Abutment

-Comes with screw

Multi-Unit Platform	SKU
4.8	131771



Multi-Unit Abutment Analog

Multi-Unit Platform	SKU
4.8	131730



Multi-Unit Abutment DIM Analog -Comes with screw

Multi-Unit Platform	SKU
4.8	133603



Multi-Unit Abutment Ti Base -Comes with screw

Multi-Unit Platform	Size	Chimney Height	SKU
4.8	RP	6 mm	131738
4.8	RP	8 mm	131740



Angled Multi-Unit Abutment Ti Base

-Comes with screw

Multi-Unit Platform	Size	SKU
4.8	RP	131734



Angled Multi-Unit Abutment Ti Base Replacement Screw

Multi-Unit Platform	Size	SKU
4.8	RP	131736


*Not a trademark of The Argen Corporation.

Multi-Unit Abutments Compatible with ASTRA TECH OsseoSpeed TX^{®*}

AO Series



Multi-Unit Abutments *-Comes with screw*

	Size	Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
	3.5/4.0 RP	4.8	1.5 mm	Str	NE	131381
	3.5/4.0 RP	4.8	2.5 mm	Str	NE	131383
	4.5/5.0 WP	4.8	1.5 mm	Str	NE	131382
	4.5/5.0 WP	4.8	2.5 mm	Str	NE	131384

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argen Corporation.

Multi-Unit Abutments Compatible with Neodent® Grand Morse™*

NGM Series



Multi-Unit Abutments *-Comes with screw*

Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
4.8	1.5 mm	Str	NE	133402
4.8	2.5 mm	Str	NE	133403
4.8	2.5 mm	Ang17	E	133404
4.8	3.5 mm	Ang17	E	133405
4.8	3.5 mm	Ang30	E	133406
4.8	4.5 mm	Ang30	E	133407

Tools

Refer to Tools section of the catalog.













*Not a trademark of The Argen Corporation.

Multi-Unit Abutments Compatible with NobelActive®*

NC Series




Multi-Unit Abutments *-Comes with screw*

	Size	Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
	3.5 NP	4.8	1.5 mm	Str	NE	131447
	3.5 NP	4.8	2.5 mm	Str	NE	131448
	3.5 NP	4.8	2.5 mm	Ang17	E	131451
	3.5 NP	4.8	3.5 mm	Ang17	E	131452
	3.5 NP	4.8	3.5 mm	Ang30	E	131455
	3.5 NP	4.8	4.5 mm	Ang30	E	131456
	4.3/5.0RP	4.8	1.5 mm	Str	NE	131449
	4.3/5.0RP	4.8	2.5 mm	Str	NE	131450
	4.3/5.0RP	4.8	2.5 mm	Ang17	E	131453
	4.3/5.0RP	4.8	3.5 mm	Ang17	E	131454
	4.3/5.0RP	4.8	3.5 mm	Ang30	E	131457
	4.3/5.0RP	4.8	4.5 mm	Ang30	E	131458



Multi-Unit Abutment Replacement Screw

	Size	SKU
	3.5 NP	131459
	4.3/5.0RP	131460

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argon Corporation.

Multi-Unit Abutments Compatible with NobelReplace® Select™*

NR Series



Multi-Unit Abutments *-Comes with screw*

	Size	Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
	3.5 NP	4.8	1.0 mm	Str	NE	131499
	3.5 NP	4.8	2.0 mm	Str	NE	131502
	4.3/5.0RP	4.8	1.0 mm	Str	NE	131500
	4.3/5.0RP	4.8	2.0 mm	Str	NE	131503
	5.0WP	4.8	1.0 mm	Str	NE	131501
	5.0WP	4.8	2.0 mm	Str	NE	131504

Tools

Refer to Tools section of the catalog.













*Not a trademark of The Argen Corporation.

Multi-Unit Abutments Compatible with Straumann Bone Level®*

SB Series



Multi-Unit Abutments *-Comes with screw*

	Size	Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
	3.3 NC	4.8	1.5 mm	Str	NE	131529
	3.3 NC	4.8	2.5 mm	Str	NE	131531
	3.3 NC	4.8	2.5 mm	Ang17	E	131533
	3.3 NC	4.8	3.5 mm	Ang17	E	131534
	3.3 NC	4.8	3.5 mm	Ang30	E	131537
	3.3 NC	4.8	4.5 mm	Ang30	E	131538
	4.1/4.8 RC	4.8	1.5 mm	Str	NE	131530
	4.1/4.8 RC	4.8	2.5 mm	Str	NE	131532
	4.1/4.8 RC	4.8	2.5 mm	Ang17	E	131535
	4.1/4.8 RC	4.8	3.5 mm	Ang17	E	131536
	4.1/4.8 RC	4.8	3.5 mm	Ang30	E	131539
	4.1/4.8 RC	4.8	4.5 mm	Ang30	E	131540



Multi-Unit Abutment Replacement Screw

Size	SKU
Universal	131541

Tools

Refer to Tools section of the catalog.

*Not a trademark of The Argon Corporation.

Multi-Unit Abutments Compatible with Zimmer Tapered Screw-Vent®*

ZS Series



Multi-Unit Abutments – Comes with screw

	Size	Multi-Unit Platform	Gingival Height	Str/Ang	Engaging/ Non-Engaging	SKU
	3.5 NP	4.8	1.0 mm	Str	NE	131594
	3.5 NP	4.8	2.0 mm	Str	NE	131596
	3.5 NP	4.8	2.5 mm	Ang17	E	131598
	3.5 NP	4.8	3.5 mm	Ang17	E	131600
	3.5 NP	4.8	3.5 mm	Ang30	E	131602
	3.5 NP	4.8	4.5 mm	Ang30	E	131604
	4.5 RP	4.8	1.0 mm	Str	NE	131595
	4.5 RP	4.8	2.0 mm	Str	NE	131597
	4.5 RP	4.8	2.5 mm	Ang17	E	131599
	4.5 RP	4.8	3.5 mm	Ang17	E	131601
	4.5 RP	4.8	3.5 mm	Ang30	E	131603
	4.5 RP	4.8	4.5 mm	Ang30	E	131605




Multi-Unit Abutment Replacement Screw

Multi-Unit Platform	SKU
Universal	131606

Tools

Refer to Tools section of the catalog.




*Not a trademark of The Argen Corporation.








ISO 13485 and MDSAP

Quality in everything is an Argen core value. ISO 13485 is a quality standard used universally around the world. Argen earned this certification by implementing a stringent quality management system that governs our end-to-end manufacturing process. MDSAP is Medical Device Single Audit Program and is the result of several countries coming together to develop a single regulatory audit that would meet all the requirements of the countries. Ensuring quality in everything, Argen is ISO 13485 certified and MDSAP compliant.

Tools

Tool	Description	SKU
	7 mm Torque Wrench	133760
	8 mm Torque Wrench	133761
	Adapter for 7 mm Torque Wrench	133762
	Adapter for 8 mm Torque Wrench	133763
	ArgenIS Screwdrivers are interchangeable with the adapter and Torque Wrench. Made of stainless steel, 24 mm in length.	
Hex 1.20 mm	Biomet 3i Osseotite® Biomet 3i Certain® MegaGen AnyRidge®	132433
Hex 1.22 mm	Dentsply Friadent Xive® Hiossen ET/HG®	132434
Hex 1.27 mm	Astra Tech OsseoSpeed TX® Zimmer Tapered Screw-Vent® BioHorizons® Astra Tech EV® CAMLOG®	132435
Torx	Straumann® synOcta® Tissue Level Straumann® Bone Level	132436
UniGrip	Nobel Brånemark System® NobelReplace® Select™ NobelActive® Multi-unit Abutment	132437
Neodent	Neodent® Grand Morse™	133759

*Not a trademark of The Argen Corporation.

Tool	Description	SKU
Insert for angled screw channel screwdriver handle 35 mm 	For insert holder and Torque Ratchet for angulated screw channel abutment screw	127562
Insert for angled screw channel screwdriver handle 24 mm 	For insert holder and Torque Ratchet for angulated screw channel abutment screw	127561
Insert for angled screw channel screwdriver handle 20 mm 	For insert holder and Torque Ratchet for angulated screw channel abutment screw	127560
Multi-Unit Abutment Key 	Multi-unit platform: 4.8	132438
Adapter for Multi-Unit Abutment Key 	For use with latch style handles and wrenches	132440

*Not a trademark of The Argen Corporation.

ArgenIS

Implant System and Scan Bodies Compatibility

Implant Info			Indication Compatibility				Scan Body Compatibility						
System	Code	Size	Ti	Hybrid	ASC	Multi-Unit Abutment	ArgenIS **	Multi-Unit ***	NT **	3Shape	Elos	Elos IO	Dess IO
Astra Tech EV	SEV	3.0, 5.4	X	X			X		X	X		X	X
	SEV	3.6, 4.2, 4.8	X	X	X		X		X	X		X	X
Astra Tech OsseoSpeed® TX*	AO	3.0 NP		X									X
	AO	3.5/4.0 RP, 4.5/5.0 WP	X	X	X	X	X	X	X	X	X	X	X
BioHorizons®*	BIO	3.0	X				X		X	X			X
	BIO	3.5 NP, 4.5 RP, 5.7 WP	X	X			X		X	X	X	X	X
Biomet 3i Certain®**	BC	3.4 NP, 4.1 RP, 5.0 WP	X	X	X		X		X	X	X	X	X
Biomet 3i Osseotite®**	BO	3.4 NP, 5.0 WP		X	X				X	X	X	X	X
	BO	4.1 RP	X	X	X				X	X	X	X	X
CAMLOG®*	CAM	3.3, 5.0	X	X			X		X	X	X	X	X
	CAM	3.8, 4.3	X	X	X		X		X	X	X	X	X
Dentsply Friadent XiVE®**	DX	3.4 NP	X		X				X	X	X		X
	DX	3.8 RP, 4.5 WP	X		X				X	X	X	X	X
	DX	5.5	X						X	X	X	X	X
Hiossen ET/HG*	HIO	3.5, 4.0/4.5/5.0/6.0	X	X			X		X				X
MegaGen AnyRidge®**	MA	Universal	X	X	X		X						X
Neodent® Grand Morse™**	NGM	Universal	X	X	X	X	X	X					X
NobelActive®*	NC	3.5 NP, 4.3/5.0 RP	X	X	X	X	X	X	X	X	X	X	X
	NC	5.5 WP	X	X			X			X			X
Nobel Brånemark System®*	NB	3.5 NP, 4.1 RP	X	X	X				X	X	X	X	X
	NB	5.1 WP	X	X					X	X	X	X	X
NobelReplace®**	NR	3.5 NP, 4.3 RP, 5.0 WP	X	X	X	X	X	X	X	X	X	X	X
	NR	6.0		X			X		X	X			X
Straumann® Bone Level*	SB	3.3 NC, 4.1/4.8 RC	X	X	X	X	X	X	X	X	X	X	X
Straumann® synOcta®*	SS	3.5 NN		X					X	X			X
	SS	4.8 RN, 6.5 WN	X	X	X		X		X	X	X	X	X
Zimmer Tapered Screw-Vent®*	ZS	3.5 NP, 4.5 RP	X	X	X	X	X	X	X	X	X	X	X
	ZS	5.7 WP	X	X	X		X		X	X	X	X	X

* Not trademarks of The Argen Corporation. ** Sold by Argen. *** Sold by Argen, one size fits all.

Ti Bases Gingival Heights

All platforms available in engaging and non-engaging

Series	Size	Gingival Heights
Astra Tech EV* (SEV)	3.6	1.0mm, 3.0mm
	4.2	1.0mm, 3.0mm
	4.8	1.0mm, 3.0mm
Astra Tech OsseoSpeed® TX* (AO)	3.0 NP	1.0mm, 3.0mm
	3.5/4.0 RP	1.1mm, 3.0mm
	4.5/5.0 WP	0.9mm, 3.0mm
BioHorizons®* (BIO)	3.5 NP	0.8mm
	4.5 RP	0.3mm
	5.7 WP	0.3mm
Biomet 3i Certain®* (BC)	3.4 NP	0.3mm, 1.5mm
	4.1 RP	0.3mm, 1.5mm
	5.0 WP	0.3mm, 1.5mm
Biomet 3i Osseotite®* (BO)	3.4 NP	0.3mm
	4.1 RP	0.3mm
	5.0 WP	0.3mm
CAMLOG®* (CAM)	3.3	0.5mm
	3.8	0.3mm, 1.5mm
	4.3	0.3mm, 1.5mm
	5.0	0.3mm, 1.5mm
Hiossen ET/HG* (HIO)	3.5	1.0mm, 3.0mm
	4.0-6.0	1.0mm, 3.0mm
MegaGen AnyRidge®* (MA)	Universal	2.0mm
Neodent® Grand Morse™* (NGM)	Universal	0.8mm, 1.5mm
NobelActive®* (NC)	3.5 NP	1.2mm, 3.0mm
	4.3/5.0 RP	1.2mm, 3.0mm
	5.5 WP	1.5mm

Series	Size	Gingival Heights
Nobel Brånemark System®* (NB)	3.5 NP	0.3mm
	4.1 RP	0.3mm
	5.1 WP	0.3mm
NobelReplace® Select™* (NR)	3.5 NP	0.3mm
	4.3 RP	0.3mm
	5.0 WP	0.3mm
Straumann®* Bone Level (SB)	3.3 NC	1.2mm, 3.0mm
	4.1/4.8 RC	1.2mm, 3.0mm
Straumann® synOcta®* (SS)	3.5 NN	0.5mm
	4.8 RN	0.4mm
	6.5 WN	0.3mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5 NP	0.8mm, 1.5mm
	4.5 RP	0.3mm, 1.5mm
	5.7 WP	0.3mm, 1.5mm

MKT-18953-Ver4-09/24

Multi-Unit Abutment

Gingival Heights and Angles

Series	Size	Multi-Unit Platform	Angle	Gingival Heights
Astra Tech OsseoSpeed® TX* (AO)	3.5/4.0	4.8	Straight	1.5 mm, 2.5 mm
	4.5/5.0	4.8	Straight	1.5 mm, 2.5 mm
Neodent® Grand Morse™*	Universal	4.8	Straight	1.5 mm, 2.5 mm
	Universal	4.8	17°	2.5 mm, 3.5 mm
	Universal	4.8	30°	3.5 mm, 4.5 mm
NobelActive®* (NC)	3.5 NP	4.8	Straight	1.5 mm, 2.5 mm
	3.5 NP	4.8	17°	2.5 mm, 3.5 mm
	3.5 NP	4.8	30°	3.5 mm, 4.5 mm
	4.3/5.0 RP	4.8	Straight	1.5 mm, 2.5 mm
	4.3/5.0 RP	4.8	17°	2.5 mm, 3.5 mm
	4.3/5.0 RP	4.8	30°	3.5 mm, 4.5 mm
NobelReplace® Select™* (NR)	3.5 NP	4.8	Straight	1.0 mm, 2.0 mm
	4.3 RP	4.8	Straight	1.0 mm, 2.0 mm
	5.0 WP	4.8	Straight	1.0 mm, 2.0 mm
Straumann® Bone Level* (SB)	3.3 NC	4.8	Straight	1.5 mm, 2.5 mm
	3.3 NC	4.8	17°	2.5 mm, 3.5 mm
	3.3 NC	4.8	30°	3.5 mm, 4.5 mm
	4.1/4.8 RC	4.8	Straight	1.5 mm, 2.5 mm
	4.1/4.8 RC	4.8	17°	2.5 mm, 3.5 mm
	4.1/4.8 RC	4.8	30°	3.5 mm, 4.5 mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5	4.8	Straight	1.0 mm, 2.0 mm
	3.5	4.8	17°	2.5 mm, 3.5 mm
	3.5	4.8	30°	3.5 mm, 4.5 mm
	4.5	4.8	Straight	1.0 mm, 2.0 mm
	4.5	4.8	17°	2.5 mm, 3.5 mm
	4.5	4.8	30°	3.5 mm, 4.5 mm

MKT-18965-Ver4-09/24

Why **ArgenIS**?

Trust

- Extensive abutment warranty – includes abutment and implant
- Easy return policy
- Committed to doing whatever it takes to make dental labs successful

Product

- Model, abutment, and crown production in one workflow
- Options for multi-unit abutments, hybrids, retention lines, and gold hue
- Narrow scan bodies allowing for more flexible use

Manufacturing

- State-of-the-art manufacturing in the USA
- Industry-leading tolerances between the implant and abutment
- FDA-cleared proprietary process and strict adherence to quality control standards

Fast One-Stop Shop

- One day to make, one day to ship
- Fast file upload and split file management

Custom Abutment and Bar Warranty

Argen provides a limited warranty for all Argen manufactured custom abutments and specified bars. In addition, if an implant supplier does not honor their warranty due to the use of an Argen abutment or bar, Argen will replace both the abutment and the implant – certain limitations apply. Please visit argen.com for complete terms and conditions of the warranty.

Visit argen.com for complete terms of the warranty.



ArgenIS Product Catalog

The Argen Corporation
8515 Miralani Drive
San Diego, CA 92126

argen.com | (800) 255-5524



Scan for Online Catalogs

