

ArgenIS

Angled Screw Channel Abutment

With the option to place the screw access hole anywhere between 0° to 25° in a 360° radius using ArgenIS Angled Screw Channel Abutments, you can now restore almost any screw-retained case without compromising esthetics.

Available Platforms

ASTRA TECH EV™*

ASTRA TECH OsseoSpeed TX™*

Biomet 3i Certain®*

Biomet 3i OSSEOTITE®*

CAMLOG®*

Dentsply Friadent XiVE®*

Megagen AnyRidge®*

NobelActive™*

NobelReplace Select®*

Nobel Brånemark®*

Straumann Bone Level®*

Straumann SynOcta®*

Zimmer Tapered Screw-Vent®*

*Not trademarks of The Argen Corporation.


Available in both engaging and non-engaging versions.

ARGEN^{IS}
Implant Systems



Angled for maximum esthetics

Featuring an **engineered surface coating** for optimal cement retention.



▶ **Learn more at** argen.com/custom-abutments

argen.com | (800) 255-5524

MKT-18493-Ver2-07/22