

THE MOST ADVANCED ZIRCONIA TECHNOLOGY

ARGEN^Z HT⁺

TRANSLUCENCY AND STRENGTH



ArgenZ HT+ is the most translucent HT material with the highest strength and blends with ArgenZ Anterior and other ceramic materials making it perfect for combination cases.

- 45% Translucency
- 1250 MPa
- Bridges up to 14 Units
- Full Contour or Frames

ARGEN^Z ANTERIOR

Super Translucent Zirconia

COSMETIC ZIRCONIA FOR THE ESTHETIC ZONE



- 50% Translucency
- 765 MPa
- Single to 3-unit Anterior Bridges
- Single Unit Posterior Crowns

ArgenZ Anterior is an exciting alternative to lithium disilicate offering the highest translucency with optimum strength most closely mimicking the natural dentin of the tooth.



Place your order online
at argen.com

ARGEN[®]

(800) 255-5524
MKT-18088-Ver1-08/19