

“Predictable esthetics, repeatable staining, and artistic results.

ArgenZ HT+ Pre-Shaded is fantastic!”

- David Jackson, Owner
Mid-South Dental Lab



TRANSLUCENCY
45%

STRENGTH
1250MPa

SHADES
21

THICKNESSES
9

ARGEN Z HT+
PRE-SHADED



Get the shades you need, **the way you want them.**

16 VITA Classical shades | 3 Bleach shades | 2 “Light” A Shades*

- + Discs or milled units
- + Pre-Shaded – For precise shade matching
- + 21 shades + white
- + Dipped – For transitionally shaded zirconia

Order discs or milled units at argen.com #ArgenAwesome

*The L or “Light” shades are designed to match the incisal third. They have less chroma than the standard A1 and A2 shades allowing for more customization with shade liquids, modifiers, or stain. Available in pre-shaded discs only.

ARGEN®
MKT-18052-Ver1-08/19

ArgenZ HT+

WHAT CUSTOMERS ARE SAYING



NEXT GENERATION HT ZIRCONIA

"ArgenZ HT+ is the best looking zirconia we've ever seen! The new pre-shaded discs are consistent and accurate to the intended hue. Our doctors are already asking for more."

- Steve Spoon, Owner,
Miracle Dental Lab

"Wow you guys nailed it just did six comparisons to other company's esthetic zirconia and by far better!!!"

- Scott Bowman,
William Scott Dental Studio

"ArgenZ HT+ Pre-Shaded achieves the intended hue and value. It's quick and easy to stain the units to match the shade tab."

- Darrick Hartman, Owner,
Professional Dental Arts Lab

"The new ArgenZ HT+ Pre-Shaded discs are consistent to the intended hue and help simplify our workflow. They create ease of use for our technicians to go from start to finish and delivery to the doctor."

- Jeremy M. Guidinger,
Digital Production Manager,
Van Hook Dental Studio

"Argen's HT+ has the perfect combination of opacity and translucency which mimics the appearance of lithium disilicate. HT+ has allowed us to increase our in-house production by designing, milling, and most importantly, sintering both hybrid abutments and crowns under one sintering program! Argen's HT+ is a total game changer!"

- Frank Escalante,
CDT AACD,
Newport Beach Smile Designs

"Shade transition from body to incisal is much better. The HT+ has a better blend which makes it much easier to nail the incisal (less blue required)."

- Robert Milligan,
CadCrowns, CA Inc.

"Much better than the competition. Argen is blowing them away with zirconia [and] doctors are giving great feedback."

- JP Phillips,
Southern Dental Restoration

"Argen's HT+ Pre-Shaded Zirconia is stronger, more translucent and a better shade match than any of the higher strength zirconia brands I tried. The shades are right on. I do a lot less layering now and my crowns are looking better than ever. By the way the customer service is terrific!"

- Mike Sandy, Owner,
Sumner Dental Lab

Order discs or milled units at argen.com #ArgenAwesome

