


EQUIPMENT & ACCESSORIES



Products you need from a partner you trust.



argen.com | (800) 255-5524



INVEST IN YOUR FUTURE

Argen empowers you with streamlined workflows integrated with industry-leading design services, our award-winning equipment and technical support, cutting-edge scanning, milling, and 3D printing technologies, and a wide array of supporting materials and accessories.



40 Member Support Team – Our equipment and technical experts are here to answer all of your questions and help you resolve any equipment issues should they arise.



Argen Swift Credit Financing – Argen Swift Credit Financing offers convenient financing options, directly from the same source that supports your digital equipment needs.



Equipment Support Program – The Argen (ESP) Equipment Support Program grants you access to over 40 equipment and technical support experts, ensuring comprehensive coverage for your software and hardware requirements.



Design Resources – Eliminate double entries and file transfer time with Argen's seamless integration with 3Shape, Evident, and EasyDent.



TABLE OF CONTENTS

Scan for Catalogs



Scanners & Accessories

3Shape Scanners.....	4
PC Options.....	8
Zero Risk Scan Box.....	10
Cintiq Monitors.....	11
Dental System.....	12
LabCare and Warranty.....	13
Accessories.....	14
Articulator Transfer Plates.....	16
Medit Scanners.....	18
Parts.....	21
exocad Software.....	22

Mills & Accessories


Roland Mills.....	24
Parts.....	26
vhf Mills.....	28
Parts.....	32
Sintering Furnaces and Parts.....	34
Suction Units.....	38
Parts.....	39
MillBox CAM Software.....	40
vhf DentalCAM Software.....	41

Printers & Accessories

Ackuretta SOL Printers.....	42
Accessories.....	42
Asiga Printers.....	44
Accessories.....	47

Argen Training

Overview.....	48
3Shape Training.....	49
Mill Training.....	50



3SHAPE SCANNERS

It starts with the scan.

Get powerful scanning ability and exceptional accuracy with 3Shape F and E Series scanners .

Design software done right.

Pair 3Shape scanners with 3Shape's best-in-class design software to create the ultimate CAD solution for your lab.

Versatile manufacturing options.

Connect to your own manufacturing equipment or send files to Argen's digital center.

Scanners & Accessories

3Shape E2 Scanner

The 3Shape E2 Scanner enhances productivity with texture scanning. The E2 scanner is compatible with the 3Shape Dental System software of your preference, allowing you to tailor it to meet your specific business requirements effectively.



CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
2x5 MP	10 μ m	<ul style="list-style-type: none">• 24 sec (full arch)• 72 sec (full arch impression)	B/W	Standard



3Shape E3 Scanner

The 3Shape E3 Scanner excels in accuracy and performance, particularly for implant bars, establishing itself as an optimal choice for dental laboratories. Driven by 3Shape technology, it delivers detail and reliability, ensuring a blend of precision and efficiency.

CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
2x5 MP	7 μ m	<ul style="list-style-type: none">• 18 sec (full arch)• 64 sec (full arch impression)	Color	Standard

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





3Shape E4 Scanner

Designed with the efficiency of dental laboratories in focus, the 3Shape E4 model can swiftly scan a full arch in a mere 11 seconds, boasting an impressive accuracy of 4 microns. This model delivers both speed and high precision, catering to the demands of high-production environments.

CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
4x5 MP	4 µm	<ul style="list-style-type: none"> • 9 sec (full arch) • 45 sec (full arch impression) 	Color	Die-in-model

3Shape F8 Scanner

Gain the power of dual model scanning with 3Shape's newest desktop scanner, the F8 is in a class of its own. Equipped with intuitive magnetic place-and-scan design, auto-start, in-software and visual lights guidance, for improved user experience.



CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
4x5 MP	4 µm	<ul style="list-style-type: none"> • Dual arch scanning • 9 sec (full arch) • 45 sec (full arch impression) 	Color	Die-in-model

Scanners & Accessories

Already have a 3Shape Scanner?

Transfer your dongle and get
industry-leading technical support.

- Live technical support
- Advanced networking expertise
- Design strategies
- Access to Argen renowned 3Shape training



The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



3Shape Scanner PC Options

High Performance PC

Lenovo P360

- Intel 12th Gen I7-12400
- Nvidia T1000 8GB
- 1TB SSD
- 32GB RAM



Ultra Performance PC

Lenovo P360

- Intel 12th Gen-12500
- Nvidia RTX A2000 6GB
- 1TB SSD
- 64GB RAM



Note: monitors are not included.

Scanners & Accessories



The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



Zero Risk Scan Box

Now you can go digital quickly and easily – no design experience or software needed. Argen’s exclusive Zero Risk Scan Box program allows you to choose from 3Shape E2, E3, and E4 Scanners*, giving you the flexibility to scan models in-house, send your cases to your preferred design service, and send to Argen for manufacturing. There is no annual subscription, and you have the option to return the scanner at any time.**



“ Argen’s Zero Risk Scan Box is an innovative solution that addresses the needs of those dental technicians who are still not comfortable with designing digitally themselves. ”

Rune Fisker, 3Shape VP Product Strategy, Ph.D., MBA, Business Unit Director, Dental Lab and Dental Clinic

- ✓ Zero Money Down
- ✓ Zero Dongle Renewal Fee
- ✓ Zero First Month’s Payment
- ✓ Zero Percent Financing
- ✓ Zero Shipping Costs
- ✓ **Zero Risk!**

Contact your Argen Sales Representative to learn more.

*36 month financing for qualified customers based on financing approval
**No refunds on payments

Scanners & Accessories

Lose the mouse and Design with Cintiq

Compatible with 3Shape and exocad CAD design solutions, unleash your creativity with your choice of Cintiq interactive pen and touch displays.

- Pressure sensitive pen
- Multi-touch glass display option*
- High definition and superb color accuracy
- Compatible with 3Shape and exocad software

* Multi touch function on applicable models



CINTIQ MONITOR	SCREEN SIZE	RESOLUTION	COLORS	MULTI-TOUCH
Cintiq 22	21.5 in (54.6 cm)	Full HD (1920 x 1080)	16.7 million	No
Cintiq Pro 22	21.5 in (54.6 cm)	Ultra HD (3840 x 2160)	1.07 billion	Yes
Cintiq Pro 24	23.6 in (60 cm)	Ultra HD (3840 x 2160)	1.07 billion colors	Yes
Cintiq Pro 24 Touch	23.6 in (60 cm)	Ultra HD (3840 x 2160)	1.07 billion colors	Yes

The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



MKT-18062-Ver3-01/24

Dental System

3Shape's Dental System is a comprehensive CAD solution that empowers you to concentrate on craftsmanship and creativity, leaving the software to handle the rest. Regardless of your lab's size, Dental System enhances productivity and facilitates the delivery of profitable services, including traditional crown and bridge, implant-retained restorations, smile design, surgical guides, clear aligners, splints, removables, and a variety of other offerings.



Complete Restorative

- Includes Dental System Premium and Dental System Removables plus all available restorative add-on modules
- One low bundled price



Premium

- Includes Dental System Crown & Bridge plus splints, digital temps, digital smile design, RealView, telescopes, post & core, and more



Crown & Bridge

- Design full anatomies, bridges, copings, frameworks, and wax-ups
- CAD points available for advanced indications



Removables

- Design full and partial dentures



Scan Only

- Receive TRIOS cases or scan gypsum models and send to Argen Design Services



Add-Ons / CAD Points

Supplement your solution with add-on modules to offer abutments, implant bars and bridges, models, dentures and more. You can also purchase CAD Points to use on a case-by-case basis for indications not included in your package.

Scanners & Accessories

Protect Your Investment

3Shape LabCare makes Dental System a future-secure investment

3Shape's LabCare subscription package provides high value services such as unlimited upgrades, support, and access to our live and recorded free educational webinars. With a LabCare subscription, your 3Shape scanner is covered by a five-year warranty and replacement services to secure your investment.









The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524










MKT-18062-Ver3-01/24

13

3Shape Scanner Accessories

ACCESSORIES	SKU	DESCRIPTION
	121907	<ul style="list-style-type: none"> • Allows the user to mount all 3Shape compatible articulators in all E-scanner models. • Requires Dental System 17.3.0 or newer Software version.
	121936	<ul style="list-style-type: none"> • Use as replacement or additional 3Shape Rubber Bands included in the delivery of the 3Shape Articulator Holder.
	115949	<ul style="list-style-type: none"> • Use for scanning of non-mounted models with all 3Shape lab scanners. • Use together with Scan Height Adapter. • For large interior scanner interface plates, please use Generic Articulator Transfer Plates.
	118874	<ul style="list-style-type: none"> • Use with D2000 scanners (left axis). • Use with interface plate on top to ensure correct position of models and dies for scanning.
	118875	<ul style="list-style-type: none"> • Use with D2000 scanners (right axis). • Use with interface plate on top to ensure correct position of models and dies for scanning.
	120018	<ul style="list-style-type: none"> • Use with and included with all E Scanners (E1, E2, E3, E4). • Use with interface plate on top to ensure correct position of models and dies for scanning.

Scanners & Accessories

ACCESSORIES	SKU	DESCRIPTION
	118876	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing the occlusal alignment of upper and lower cast during scanning.
	115828	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing complete wax-up bridges in the scanner.
	114865	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing dental impressions in the scanner.
	118878	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing implant models in the scanner.
	118879	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing full arch single and triple tray impressions in the scanner.
	122224	<ul style="list-style-type: none"> • Fixture used for fixing up to 7 individual dies. • Use together with scan height adapter with large interior lab scanners.
	115822	<ul style="list-style-type: none"> • 3Shape Articulator used for articulation of digital lab models designed in 3Shape's Model Builder.








The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524









3Shape Articulator Transfer Plates

You need both an Articulator Calibration Object and the Articulator Transfer Plates to transfer the model's position from physical articulator to virtual articulator in Dental System.

ACCESSORIES	SKU	DESCRIPTION
 <p>Articulator Calibration Object</p>	118882	<ul style="list-style-type: none"> • Fixture for all desktop scanner models. • Use for fixing implant models in the scanner.
 <p>Generic Articulator Transfer Plate - 1 stripe</p>	120298	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis). • Allows users to make their own transfer plates for articulator systems – not (yet) supported by other 3Shape Transfer Plates, such as BioArt & many more.
 <p>Generic Articulator Transfer Plate - 2 stripe</p>	120299	<ul style="list-style-type: none"> • Use with D2000 scanners (right axis). • Allows users to make their own transfer plates for articulator systems – not (yet) supported by other 3Shape Transfer Plates, such as BioArt & many more. • Use in conjunction with Generic Transfer Plate – 1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential.
 <p>Adesso Artex Articulator Transfer Plate - 1 stripe</p>	118883	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis). • Use with Baumann Adesso coupling system typically used with Ivoclar Stratos 300. • Also used with ACR (Artex compatible) articulators.
 <p>Adesso Artex Articulator Transfer Plate - 2 stripe</p>	118884	<ul style="list-style-type: none"> • Use with D2000 scanners (right axis). • Purchased in addition to Adesso Artex Articulator Transfer plate – 1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential. • Use with Baumann Adesso coupling system typically used with Ivoclar Stratos 300. • Use with ACR (Artex compatible) articulators.
 <p>Whip Mix Articulator Transfer Plate - 1 stripe</p>	118885	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis). • Use with Whip Mix DENAR Mark 330 articulators.
 <p>Whip Mix Articulator Transfer Plate - 2 stripe</p>	118886	<ul style="list-style-type: none"> • Use with D2000 scanners (right axis). • Use in conjunction with Whip Mix DENAR Mark 330 Transferplate –1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential.

Scanners & Accessories

ACCESSORIES	SKU	DESCRIPTION
	118888	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis)). • Use with SAM articulators.
	118889	<ul style="list-style-type: none"> • Use with SAM articulators and D2000 scanners (right axis). • Use in conjunction with SAM Transfer plate – 1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential.
	118891	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis)). • Use with KaVo Protar articulators.
	118892	<ul style="list-style-type: none"> • Use with KaVo Protar articulators and D2000 scanners (right axis). • Use in conjunction with KaVo Protar Transfer plate – 1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential (see note 2).
	124403	<ul style="list-style-type: none"> • Use with all large interior scanners (E1, E2, E3, E4, and D2000 (left axis)). • Use with Shofu articulators. • Requires MSMP plate to be mounted on top. Please purchase MSMP plate from your articulator provider.
	124404	<ul style="list-style-type: none"> • Use with Shofu articulators and D2000 scanners (right axis). • Use in conjunction with Shofu Transfer plate – 1 stripe (left axis) to be able to benefit from the D2000 scanner's full potential. • Requires MSMP plate to be mounted on top. Please purchase MSMP plate from your articulator provider.

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





T-Series lab scanners speed up your workflow.


CAD/CAM dental work requires the highest accuracy to produce well-designed bridges, implants, and bars. The Medit T-Series ensures accurate, high-quality scans that adhere to your laboratory's standards. The T-Series adapts to your workflow with efficiency, by using the flexible multi-die to scan a full-arch or partials with multiple dies simultaneously. Auto-elevation, wider scan area and no blind spots are some of the highlights of the T-Series.

Customized scan sequence.

Enjoy the flexibility of choosing your own scan sequence with Medit Scan for Labs. Simply drag and drop the scan stages and reorganize the sequence to suit your working style and needs!


Scanners & Accessories

Medit T-Series




Super Fast Scanning

With our high-quality hardware and optimized software, you can scan a full-arch in just 8 seconds with the T710, Medit's fastest dental lab scanner.



High-resolution Cameras

5.0MP cameras ensure high-resolution detailed scan data. With the four-camera system, the T710 lab scanner covers a wide scan area, eliminating any blind spots.



Auto-elevation

The auto-elevation feature allows you to say goodbye to stacking jigs to get the correct scanning position.



High Accuracy

4-micron accuracy



Open System

Enjoy the freedom of Medit's open system which allows you to import and export files in STL format so you can design on virtually any software.



Medit 710 Scanner

From the company that brought you the first blue-light tabletop scanner, introducing the Medit T710, the fastest Medit tabletop scanner you have yet to experience.

CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
4x5 MP	4 μ m	<ul style="list-style-type: none"> 8 sec (7 cut) 45 sec full arch 	Color	Phase-shifting optical triangulation

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



Medit 510 Scanner

The advanced, high speed positioning system of the new T Series is designed for optimal performance for your laboratory, speeding up your workflow and increasing productivity.



CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
2x5 MP	7 μ m	• 12 sec (7 cut)	Color	Phase-shifting optical triangulation



Medit 310 Scanner

Scan more objects at the same time thanks to the wider scan area of our T-Series tabletop scanners. Capture vivid color texture scans, including hand-marked margins or handwritten notes.

CAMERAS	ACCURACY (ISO 12836)	SCAN SPEED	TEXTURE	SCANNING STRATEGY
2x5 MP	9 μ m	• 18 sec (7 cut)	Color	Phase-shifting optical triangulation

Scanners & Accessories

Medit Scanner Parts

<u>PART</u>	<u>SKU</u>
Calibration Panel	133321
LED Protector	133322
Multi-Die Jig	133323
Single Die (16x)	133324
Spray Supporter	133325
Screw Jig	133326
Medical Power Adapter	133327
Power Cable	133328
USB Cable	133329
KAS Jig	133330
AM Jig	133331
Articulator Plate	133332
Impression Scanning Jig	133333

The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524





exocad SOFTWARE

Powerful CAD solution for dental labs.

exocad CAD software offers swift operation and a user friendly design, which reduces training expenses while enhancing productivity. It handles intricate cases reliably and efficiently. exocad's DentalCAD® is designed using cutting-edge research on organic modeling, ensuring speed and flexibility. exocad's dedicated engineers design their software based on customer feedback, keeping exocad users ahead of industry trends.


















Extensive, versatile, integrated.

Effortlessly expand your services with exocad's integrated add-on modules in DentalCAD. Enjoy a smooth, stable workflow, increased functionality, and cost-efficiency tailored to your needs.

Scanners & Accessories

A solution that grows with your business: Select between standalone DentalCAD modules or powerful Bundles, as Perpetual or FLEX license options.

exocad Packages for DentalCAD

	Product	Core Lab Version	Advanced Lab Bundle	Implant Lab Bundle	Ultimate Lab Bundle
	DentalCAD Core Version	✓	✓	✓	✓
	Virtual Articulator	Add-on	✓	✓	✓
	Provisional Module	Add-on	✓	✓	✓
	TruSmile Module	Add-on	✓	✓	✓
	ZRS Tooth Library	Add-on	✓	✓	✓
	Implant Module	Add-on	Add-on	✓	✓
	Bar Module	Add-on	Add-on	✓	✓
	DICOM Viewer	Add-on	Add-on	✓	✓
	Model Creator	Add-on	Add-on	Add-on	✓
	Smile Creator	Add-on	Add-on	Add-on	✓
	FullDenture Module	Add-on	Add-on	Add-on	✓
	PartialCAD Module	Add-on	Add-on	Add-on	✓
	Bite Splint Module	Add-on	Add-on	Add-on	✓
	Jaw Motion Import	Add-on	Add-on	Add-on	✓
	xSNAP Module	Add-on	Add-on	Add-on	Add-on
	In-CAD Nesting	Add-on	Add-on	Add-on	Add-on
	Nesting (includes In-CAD Nesting)	Add-on	Add-on	Add-on	Add-on

exocad MKTPL-DentalCAD-2023-02-en

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





Roland Mills

Milling made easy.

Get ultimate flexibility and productivity thanks to Roland's sophisticated technology and ability to mill a variety of materials.

Your mill. Your way.

Mix and match your preferred mill, suction unit, furnace, and accessories to meet your specific needs.

Your all-access pass.

Purchase a mill through Argen and get access to our tooling, milling strategies, and renowned technical support team.

Mills & Accessories




Roland Mill DWX-52D Plus

The Roland DWX-52D Plus dry mill is the next generation of the industry's best-selling dental mill, making the mill you know and love even better. The DWX-52D Plus offers better disc handling, more material capabilities, improved tool management and other enhanced features, to deliver even more efficiency.

MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Wax • PMMA • CoCr 	1	15	Brushless DC Motor	6,000 to 30,000 rpm

Roland Mill DWX-52DCi

The Roland DWX-52DCi dry mill is the ultimate in ease, material flexibility and unattended production. With a 6-slot Automatic Disc Changer, 15-station Automatic Tool Changer and other automated features, it gives users of all skill levels an exceptionally powerful and efficient dental restoration milling solution.



MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Wax • PMMA • CoCr 	6	15	Brushless DC Motor	6,000 to 30,000 rpm

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





Roland Mill DWX-53DC

The DWX-53DC revolutionizes dental labs with its unmatched technology: a redesigned spindle with 3x gripping power, up to 4mm tooling for rapid PMMA roughing, and exclusive features. It includes a built-in webcam for real-time monitoring, either on-site or remotely, and comes with a dedicated extraction unit.

MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Wax • PMMA • CoCr 	6	15	Brushless DC Motor	6,000 to 30,000 rpm

Roland Mill Parts

ACCESSORIES

	SKU
DWX-50 Clamp Gasket	129439
DWX-50 Pin Type Mat Adap	124341
DWX-50 Bur Holder	124329
DWX-50 Tool Setter	123229
DWX-50 Disc Clamp	119708
DWX-51D Clamp Screw	121050
DWX Collet	123330
DWX-52 Disc Adapter (6)	122140
DWX-52 Disc Adapter Rack	122169
DWX-52D Time Reduc Kit - Model	124297

ACCESSORIES

	SKU
DWX52D/DCi Pin Mat Adapt	124346
DWX TRK Denture	124429
DWX TRK Denture Tips	124430
DWX TRK Denture Jig	124431
Adapter Barcode A-F	124514
Adapter Barcode G-Z	124515
DWX-52D C Clamp	129433
DWX-53 Disc Adapter (6)	132767
DWX-53DC C Clamp	132768
DWX-53 Disc Adapter Rack	132769

Mills & Accessories

Roland Mill Parts

ACCESSORIES

	SKU
DWX-53DC Collet	132773
Coolant 1000ml	119024
DWX-4W Jaeger Collet	124345
Chelating Agent (250ml)	127021
DWX-42W Collet	124432
DWX-42W Collet Brush	127024
DWX-42W Grease	127025
DWX-42W Collet Cleaner	127026
DWX-42W Nozzle Brush10pk	127027
DWX-42W Tool Remover Jig	127028
Coolant Filter Gen3 (2pk)	129393
CoolantFilterCapGen2 2pk	129394

REPAIR PARTS

	SKU
DWX-50 Spindle	118841
DWX-50 Detection Pin	119707
DWX-50 ATC Motor	120819
DWX-50 Air Regulator	120954
DWX-50 Bearing DU-24H	121300
DWX-50 Pulley DU-22H	121301
DWX-50 Z Wire Gear Drive	121746
DWX-50 Z Limit Switch	122246
DWX-50 Sensing Spring	123327
DWX-51D Spindle Belt	119884
DWX-4 Spindle Kit	119885
DWX-4 ATC	124343
DWX-51D X/Y Axis Screw	124504

REPAIR PARTS

	SKU
DWX-51D Lower Base Clamp	125265
DWX-52D Clamp Lever 1	129601
DWX-52D Clamp Lever 2	129602
DWX-52D Base Clamp	131747
DWX-52DCi Tray Plate	128894
DWX-52DCiChngrDoorSwitch	132660
DWX-53DC Spindle Belt	132771
DWX-53DC Spindle Kit	132772
DWX-4W Sensing Pin	121891
DWX-42W Door Switches	133601
DWX-42W Outer CheckValve	128158
DWX-42W Coolant Cap	129438

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





vhf Mills

Maximum precision.

Work with everything from wax, plastics, PEEK, zirconia and composites to cobalt-chrome-based NEM materials and even titanium and glass ceramics – quickly and with utmost accuracy.

Durable and precise.

vhf tool specialists continually engage in research and development in collaboration with material manufacturers and laboratory users. This collaborative effort aims to produce the most effective tools possible.

Efficient and innovative.

Save valuable time with vhf's technical innovations that allow you to work even more efficiently. No matter what you have in mind, the E, K, S and N Series mills support you optimally through the machining process, producing accurate dental restorations and appliances.

Mills & Accessories



vhf E4

With a machine weight of only 61 lbs and no use of compressed air, the E4 offers maximum flexibility in terms of installation. You can set the E4 anywhere. The E4 makes the milling of glass ceramics, zirconia, composites, PMMA, and wax a reality with an easy-to-implement workflow. The vhf E4 is true to the motto: unpack, connect, start milling.

MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Glass Ceramics • Zirconia • Composites • PMMA • Wax 	1	7	High Frequency 800 Watt Spindle	60,000 RPM

vhf E5

The E5, powered by vhf's innovative AIRTOOL technology, eliminates the need for compressed air. This feature provides maximum flexibility in selecting the installation site and ensures minimal operating costs. The E5's open system architecture facilitates a swift and easy transition into digital dental restoration production by seamlessly integrating into your workflows. With integrated CAM software, you can commence operations without delay.



MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Composites • PMMA • Wax • CoCr Sintering Metals 	1	16 + 1 AIRTOOL	High Frequency 800 Watt Spindle	60,000 RPM

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524





vhf K5+

The new K5+ allows you to achieve perfect milling results. This mill processes a wide range of materials and indications (blanks up to 40 mm) with five axis. It also has revolutionary technologies, such as tool-free blank clamping and an integrated ionizer, which help you achieve high-end results. The particularly stable machine bed construction made of a solid cast body ensures first-class machining results.

MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Wax • PMMA • CoCr 	1	16	High Frequency 500 Watt Spindle	60,000 RPM

vhf S5

The S5 is a highly versatile dental milling machine. It has five simultaneously operating axes, a blank changer for eight blanks and is designed for both dry and wet machining. With the optional wet grinding module, you can grind all common types of glass ceramic materials with the same high-volume production experienced in dry-milling mode.



MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Zirconia • Wax • PMMA • CoCr • Titanium • Glass Ceramic 	8 Disc Adapters Up to 24 Blocks Up to 48 Abutments	16	High Frequency 600 Watt Spindle	60,000 RPM

Mills & Accessories



vhf N4+

Designed for the fabrication of crowns, bridges, inlays, onlays, abutments, and veneers, the N4+ is an ideal choice for dental laboratories aiming for efficient, high-quality results. Combining the N4+'s wet grinding capabilities with a dry milling machine allows users to operate synergistically, optimizing their laboratory's capabilities. The N4+ ensures outstanding precision, resulting in Ultra HD restorations.

MATERIALS	HOUSED ADAPTERS	HOUSED BURS	SPINDLE MOTOR	SPINDLE SPEED
<ul style="list-style-type: none"> • Glass Ceramics • Zirconia • Composites • Plastics • Titanium • CoCr 	3	8	High Frequency 800 Watt Spindle	80,000 RPM

vhf Mill Parts

ACCESSORIES	SKU
vhf PowerSupply Cord USA	131074
KE/SE Insert Toolchanger	131077
KE/SE Calibration Blank	131081
KE/SE Collet Chuck 3mm	131086
vhf Spindle Ser Set	131087
vhf Drill 3mm Shank	131091
vhf KE+ Ivotion Kit	133226
SE USB Cable	131008
SE Blank Frame	131020
SE Calibration Set	131050
SE Maint Kit S Model SLO	131059
SE Torque Wrench 1.8Nm	131070
SE Unlocking Tool	131071
vhf Wet Grinding Opt NO1	132540
vhf Inter Brush 0.45mm	132545
vhf Inter Brush 0.6mm	132546
vhf FineFilter Cartridge	132548
vhf Flat Fold Filter	132549
vhf PH Test Strips	132550
vhf Measuring Cup 500ml	132551
vhf Coarse Filter Mat	132552
vhf Tec Liquid Pro 1L	132553
vhf Tec Powder (10PU)	132554
vhf SE Self Lock Frame	132555
vhf SE BlankFrameTwr w/8	132557
vhf SE Block Holder	132558
vhf Dess Abutment Holder	132735
vhf NTTrading AbutHolder	132736

ACCESSORIES	SKU
vhf Medentika AbutHolder	132737
vhf E5 Ivotion Kit	133304
vhf E5 Block Holder	133305
vhf E5 Bur Starter Kit	133309
vhf E5 Blank Holder	133486
vhf E5 FixingTool Insert	133487
vhf E5 AirTool Insert	133488
vhf E5 Tool Mag Insert	133489
vhf E5 Tool Mag Rubber	133493
vhf E5 Maintenance Kit	133495
vhf E5 Spindle Collet	133498

REPAIR PARTS	SKU
KE+ Glass Fuse T8 (10pc)	130965
KE+ Bellows	130973
KE+ Folding Bellow X	130976
KE+ Folding Bellow Y	130977
KE+ Nozzle Cover	130980
vhf Glass FuseT6 (10pc)	131073
KE/SE Hinged Door Lock	131076
SE Housing Foot	131014
SE Blank Frame	131020
KE/SE Hinged Door Lock	131076
SE Maint Kit S Model SLO	131059
SE Sealing Tape	131060
SE Measuring Pin	131068
vhf Membrane Pump	132543
SE Unlocking Tool	131071

Mills & Accessories

vhf Mill Parts

REPAIR PARTS

REPAIR PARTS	SKU
vhf E5 StartButton Board	133463
vhf E5 Connector Set	133464
vhf E5 Hall Sensor	133465
vhf E5 IO Board w/Cables	133466
vhf E5 LED Chamber Light	133467
vhf E5 DoorLock Solenoid	133468
vhf E5 Measuring Key	133469
vhf E5 Spindle Cable	133470
vhf E5 Control Unit	133471
vhf E5 Cover Latch	133472
vhf E5 Door Gasket	133473
vhf E5 Housing	133474
vhf E5 Lower Cover	133475
vhf E5 Protective Cover	133476
vhf E5 Upper Cover	133477
vhf E5 Door Stop	133478
vhf E5 Chamber Door	133479
vhf E5 A/B Axis	133480
vhf E5 X Axis Ball Screw	133481
vhf E5 Y Axis Ball Screw	133482
vhf E5 Z Axis Ball Screw	133483
vhf E5 Spindle Bellows	133484
vhf E5 Bellows Coupler	133485
vhf E5 Locking Bolt	133490
vhf E5 Key Protect Cap	133491
vhf E5 Tension Spring	133492
vhf E5 Emergency Release	133494
vhf E5 Packaging Set	133496

REPAIR PARTS

REPAIR PARTS	SKU
vhf E5 Transport Lock	133497
vhf E5 Spindle	133499
vhf E5 OpenEnd Wrench	133500
vhf E5 Torque 20 Bit	133501

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



Sintering Furnaces


With over 30 years of experience, Shenpaz incorporates their industry knowledge to design and manufacture premium dental furnaces with industry-leading technology. Shenpaz furnaces offer unlimited heating and cooling stages, with holding duration choices to provide various programming capabilities. These furnaces are universally compatible with a wide variety of dental zirconia

SINTRA PRO 90zr and 120zrf are supplied with one sintering tray.

Shenpaz Sintra Pro 90zr

SINTRA PRO 90zr capacity from 30 units per tray to 90 units in 3 stackable trays.


- 1 furnace unit
- 1 power cable
- 1 firing table
- 1 sintering crucible
- 1 sintering crucible top
- 1 pair of tweezers
- 1 USB: user manual
- 1 lift key



Shenpaz Sintra Pro 120zrf

SINTRA PRO 120zrf capacity from 40 units per tray to 120 units in 3 stackable trays and full discs of 98mm.

- 1 furnace unit
- 1 power cable
- 1 firing table
- 1 sintering crucible
- 1 sintering crucible top
- 1 pair of tweezers
- 1 USB: user manual
- 1 lift key



Mills & Accessories

Shenpaz Sintering Furnace Parts

<u>PARTS</u>	<u>SKU</u>
Sintra RM Beadless Tray	121500
Sintra Pro 90 Tray Set - 1 Set	127748
Sintra Pro 120 Tray Set - 1 Set	127751
Sintra Plus Fuse - 20 Amp	118907
Sintra Plus Fuse - 25 Amp	119928
Sintra Thermocouple	121294
Sintra Plus lift Cone	121457
Sintra Cross Flow Fan	123323
Sintra Heating Element - Individual	124339
Sintra Ceramic Holder	125234
Sintra Contact Strap	125235
Sintra Plus Lift Motor	127869
Sintra Pro 90 Heating Element Kit - 3 Pack	127749
Sintra Pro 120 Heating Element Kit - 4 Pack	127753
Sintra Pro 90 Lift Cone	127795
Sintra Pro 120 Modular Table	127750
Sintra Pro 120 Table Insert	127752
Sintra Calibration Disc - Individual	122225
Sintra Calibration Disc - 15 Pack	132652
Sintra Beads - 100g	124340

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



KDF Zircom

KDF Zircom Speed

Zircom and Zircom Speed sintering furnaces enable the storage of 100+ programs of 20 steps each, which cater to a diverse range of zirconia applications. Equipped with a molybdenum disilicide heater, these models ensure clean, dependable, and uniform results with every use.

The Zircom Speed, equipped with software and specialty trays allows for fast sintering cycles.

- Setter for high-speed sintering (Speed only)
- Firing plate (Speed only)
- Quick firing tray
- Firing plate (speed only)
- Quick firing stand
- Quick firing tray stand
- Quick firing cover
- Firing plate (speed only)
- Firing beads



KDF Zircom Plus

The Zircom Plus sintering furnace offers 100+ programs of 20 steps each for a wide range of zirconia. It also features a molybdenum disilicide heater that provides clean, reliable, and consistent results

- Quick firing stand
- Quick firing tray
- Quick firing tray stand
- Quick firing tray stand
- Quick firing cover
- Firing beads



Mills & Accessories

KDF Sintering Furnace Parts

SINTERING FURNACE PART	SKU
Zircon Firing Stand	130542
Zircon Firing Tray Stand	130541
Zircon Firing Trays - Slotted	129391
Zircon Firing Tray Cover	130543
Zircon High Speed - Firing Setter	130548
Zircon High Speed - Firing Tray	130547
Zircon High Speed - Firing Plate	130546
Zircon Plus Lower Firing Stand	131275
Zircon Plus upper Firing Stand	131274
Zircon Plus Firing Tray Saucer	131276
Zircon Plus Firing Tray - Slotted	131277
Zircon Plus Firing Tray Cover	131278
Zircon Thermocouple	130407
Zircon Plus Thermocouple	131272
Zircon Beads	130544
Zircon Heating Element - Individual	130545
Zircon Heater PCB	132651
Zircon Plus Heater PCB	131273
Zircon Lift Motor Assembly	132366
Zircon Fan Filter - 5 Pack	132464

The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



Suction Units

Quatro Power Station

- Quatro Power Station 2 – Brushed
- Quatro Power Station 2 – Brushless

The Quatro Suction Units are automatic dust collectors with high capacity, sealable, and easily disposable filter bags. Both models provide users with a small operation footprint.

- 220 CFM
- Brushed and Brushless Options



Quatro iVac Twin

- Quatro iVac Twin – Brushless

Designed for systems that generate large volumes of dust, the iVAC is an intelligent “next generation” dust collector with FCS (Filter Cleaning System).

- 282 CFM
- Brushless

BOFA DustPRO

BOFA’s DustPRO Universal unit filters swarf and particulates during the milling process. Its durable, brushless high vacuum pump runs continuously for extended periods without the need for frequent servicing.

- 167 CFU
- Brushless
- Dual Filtration



Mills & Accessories

Suction Unit Parts

QUATRO PARTS & ACCESSORIES

	SKU
PS2 Main Filter Bag - 2 Pack	130542
PS3 Main Filter Bag - 2 Pack	130541
PS6 Main Filter Bag - 2 Pack	129391
Drum Filter iVac Twin - Legacy (Pre 2019)	130543
Drum Filter iVac Twin - Current	130548
PS2 HEPA Filter	130547
PS3 HEPA Filter	130546
PS6 HEPA Filter - Requires 2	131275
iVac Twin HEPA Filter	131274
AR175 Brushed Motor	131276
iVac Shaker Motor	131277
iVac Isolator Kit	131278

BOFA PARTS & ACCESSORIES

	SKU
Main Filter Bag	118828
HEPA Filter	118829
Replacement Motor	125333



The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



MillBox CAM Software

MillBox is the dental CAM solution developed for milling any kind of restoration with any material. Equipped with a easy-to-follow user interface that is both innovative and appealing, MillBox simplifies the nesting of dental restorations for efficient toolpath creation.

- Import any type of dental restoration from any open CAD source.
- Automatic features are applied to each imported object: support pins, margin line detection, offset, and more.

ARGEN EDITION	SKU
Argen Roland Edition	132804
Argen vhf Edition	132803

MILLBOX ADD-ONS	ITEM
Advanced Module	125261
Strategy Editor	132805
Add vhf	132429
Add Roland	132430



Mills & Accessories


vhf DentalCAM Software

Software development has been an integral part of vhf's DNA for over 30 years. Perfectly matched to vhf milling machines, DentalCAM software offers you an intuitive way of working thanks to its clear user guidance and consistent three-dimensional presentation which help you achieve excellent results.

- Open STL import feature and utilizes optimized CAM strategies for the fastest possible processing.
- All relevant data on blanks, clamping fixtures and tools are stored in the CAM software and are regularly updated.
- Matched to efficient and high-precision tooling, designed and manufactured by vhf.

vhf DENTALCAM	SKU
DentalCAM 8 N Series	133313
DentalCAM 8 K Series	132730
DentalCAM 8 S Series	132731

vhf DENTALCAM UPGRADES	SKU
DentalCAM 6 to DentalCAM 8	132732
DentalCAM 7 to DentalCAM 8	132733



The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



Ackuretta SOL Printers

SOL is an LCD 3D printer with a monochromatic print screen and is specifically designed for same-day results in your laboratory. It prints at three times the speed of SLA printers and has five times more UV Power than earlier models so that it can match DLP Printers. SOL is your versatile solution for modern digital dentistry.

ACKURETTA SOL

- Technology: LCD Monochrome
- Wavelength: 385 – 405nm
- Accuracy: 49µm
- Build Plate Size: 128mm (5.03") x 80mm (3.15")
- Speed 45mm per hour
- Open Material System
- 2 Year Warranty*

**Must purchase Ackuretta Concierge Services*



Ackuretta SOL Printer Accessories

<u>ACKURETTA PRINTER PART</u>	<u>SKU</u>
<u>SOL Build Plate - Large</u>	<u>132565</u>
<u>SOL Build Plate - Medium</u>	<u>132707</u>
<u>SOL Build Plate - Small</u>	<u>132708</u>
<u>SOL Vat Tray</u>	<u>132566</u>
<u>SOL Print Screen</u>	<u>132567</u>

Printers



ACKURETTA CURIE - UV

- Curing Space 5.6 (L) x 5.7 (W) x 3.2 (H) in
- SD Card Reader
- Accuracy and Consistency – Validated
- Wavelength of LED – 365 / 385 / 405nm
- Dimensions – 7.9 (L) x 8.7 (W) x 6.3 (H) in

Accessories Included:

- 1 x Power Cord
- 1 x SD Card
- 1 x SD Card Reader
- 1 x Glass Plate
- 1 x Quick Start Guide



ACKURETTA CURIE PLUS - UV

- Curing Space 7.9 (L) x 5.9 (W) x 2.4 (H) in
- SD Card Reader
- Accuracy and Consistency – Validated
- Wavelength of LED – 365 / 385 / 405nm
- Dimensions – 14.4 (L) x 10.2 (W) x 9.2 (H) in
- N2 Compatible

Accessories Included:

- 1 x Power Cord
- 1 x SD Card
- 1 x SD Card Reader
- 1 x Glass Panel Box
- 1 x Ethernet Cable
- 1 x WiFi Dongle
- 1 x Quick Start Guide



NK OPTIK OTOFASH

- Curing Space 4.72 (L) x 4.7 (W) x 2.0 (H) in
- Wavelength of LED – 280 - 580nm
- 10 flashes per second
- Dimensions – 13 (L) x 13 (W) x 5.5 (H) in
- N2 Compatible

Accessories Included:

- 1 x Power Cord
- 1 x 280N2 Tray
- 1 x 360N2 Tray
- 1 x Filter (installed)
- 1 x Quick Start Guide



ACKURETTA: CLEANI

- Portable and Compact
- Build Platform Wash
- Dual Tanks – Dirty Bath, Clean Bath*
- Accelerated Workflow

Accessories Included:

- 1 x Power Supply and Cable
- 2 x Universal Jig
- 2 x Mesh Baskets
- 4 x Bars
- 1 x Quick Start Guide

**Does not include Isopropyl Alcohol*

The Argen Corporation

argen.com | customersuccess@argen.com | (800) 255-5524



Asiga MAX UV Printer

The Asiga MAX is an advanced digital dental 3D printer offering exceptional productivity in a small footprint. With 62µm HD print precision, the Asiga MAX is optimized for orthodontics, crown and bridge, surgical guides, dental models, custom trays, and dentures in lab and clinical environments. All Asiga printers are open to an extensive library of validated materials for maximum flexibility and economy.

Asiga MAX UV

- Lifetime technical support guarantee
- WiFi connectivity
- Automated support generation
- Touchscreen interface
- Proven manufacturing reliability



TECHNICAL DATA

Pixel size X,Y	62 µm
Build size X Y, Z	119 x 67 x 75mm; 7.83 x 2.64 x 3 inches
Light source	385nm
Z resolution	Variable in 1 µm
Material system	Open, use any 3rd party material
File inputs	STL, SLC, STM
Software	Composer included (lifetime software updates included)
Network compatibility	WiFi and ethernet
Industry sectors	Digital dentistry / lab production
System size	260 X 380 x 370mm
System weight	17.5kg
Packaged size/weight	410 X 500 x 480mm / 20kg
Power	12VDC 10a
Warranty	12 months manufacturer warranty

Printers



The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



Asiga Pro 4K Printer

The PRO 4K combines the latest 4K imaging technology with Asiga's proven Smart Positioning System (SPS™*) to produce a build volume 3x larger than desktop 3D printers. Precision, reliability, speed and an Open Material System all come standard to provide production continuity for the most demanding digital manufacturing environments.

The PRO 4K incorporates pixel shifting technology to increase part accuracy and surface definition without affecting printing time. 4K mode is a standard feature on every PRO 4K.

* Not a trademark of The Argen Corporation.



TECHNICAL DATA

Build Volume X, Y, Z	217 x 122 x 200mm; 8.54 x 4.8 x 7.87 inches
Pixel size - 4K mode	56µm
Pixel size - native mode	80µm
Technology	DLP
LED wavelength	365nm, 385nm or 405nm
Material compatibility	Open Material System. Over 400 validated materials available via Asiga's Material Library online.
Production	Dental, Audiology, Jewellery, Precision Additive Manufacturing
Software	Asiga Composer software. Lifetime updates included.
File inputs	STL, PLY, SLC, STM (STM - Asiga Stomp file format)
Network compatibility	WiFi and ethernet
Power requirements	100-240VAC, 50/60Hz, 500 Watts (100V - 5Amp Max. 240V - 2.1Amp)
System sizing	465 x 540 x 1370mm / 140 kg; 18.3 x 21.2 x 53.9 inches / 309 lb
Packed sizing	900 x 700 x 1540mm / 205 kg; 35.4 x 27.6 x 60.6 inches / 452 lb
Warranty	12 months manufacturers warranty
Technical support	Unlimited lifetime technical support included
Bundle includes	3D printer, Composer software, 1Kg Asiga material, 2L build tray, Asiga Flash post-curing chamber, calibration toolkit

Printers

Asiga Printer Accessories

ROLAND PRINTER PART	SKU
Asiga Max Build Platform	127804
Asiga Pro Build Platform	127809
Asiga Max Build Tray 1L	127806
Asiga Max Build Tray 2L	127807
Asiga Max Build Tray 5L	127808
Asiga Max Build Tray 10L	127805
Asiga Max UltraGloss Tray 1L	133436
Asiga Pro Build Tray 2L	127812
Asiga Pro Build Tray 5L	127813
Asiga Pro Build Tray 10L	127810
Asiga Pro UltraGLOSS Tray 2L	133438
Asiga Max Tray Storage	128521
Asiga Pro4K Tray Storage	133278

The Argen Corporation
 argen.com | customersuccess@argen.com | (800) 255-5524





EXPERT TRAINING

Elevate your skill set and grow your business.

The key lies in our exceptional trainers. Argen's CAD/CAM experts are renowned in their field, equipped with extensive real-world lab experience. Our devoted training team is committed to enhancing your technical skill set. Whether you aspire to operate a crown-a-week pace or establish a comprehensive digital lab, we are here to help you progress with confidence.



Argen Training

3Shape Training

Basic Training

Argen's comprehensive training program encompasses essential aspects such as system setup, maintenance, scanning, and design for both fixed restorations and removables. Tailored to meet the evolving needs of our customers, our webinars regularly feature diverse topics. From introducing innovative design concepts to guiding participants on digitally crafting traditional restorations, we ensure a dynamic and relevant learning experience for all.

BASIC TRAINING	LOCATION	ATTENDEE LIMIT	COST*
Webinar	Online	Unlimited	Free
Argen HQ Training	San Diego, CA	2 labs; 2 attendees per lab	\$1,000 per day
In-Lab Training	Lab Location	At lab's discretion	\$2,000 per day
3Shape HQ Training	Warren, NJ	Pending 3Shape availability	\$1,000

*Customer is responsible for all travel and accommodation fees of the lab or trainer, if applicable. New equipment purchase from Argen may be required to qualify for training.

The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



Mill Training

Choose from two training options based on your specific needs and learn how to make your mill work for you.

	Webinar
Duration	4 hours
Location	Online
Attendee Limit	2-4
Cost*	Free
CONTENT	
Mill set up	✓
Connecting CAM PC	✓
Connecting Suction Unit	✓
Mill Tests	✓
Nesting Software	✓
Milling Units	✓
Mill Maintenance	✓
Post Processing of Units	
Sintering Oven Setup	
Scanner Setup	
Design Software	

*Customer is responsible for all travel and accommodation fees of the lab or trainer, if applicable. New equipment purchase from Argen may be required to qualify for training.

Argen Training

“ We made the decision a few years ago to transfer our 3Shape licensing over to Argen based exclusively on the positive comments we heard about their 3Shape tech support...We scheduled several full day and half day webinars to assist us in some basic 101 training as well as some advanced training. When the training was completed, we scheduled the Argen trainer for a 3-day on-site visit to our laboratory – twice. The Argen trainer spent time with our scan and design technicians and also assisted us in fixing some set-up and software glitches. Argen’s training team continues to be a valuable asset to our company through Argen 3Shape support. ”

**Steve Edmonds,
Edmonds Dental Co., Inc**



The Argen Corporation
argen.com | customersuccess@argen.com | (800) 255-5524



Argen Equipment Product Catalog

The Argen Corporation
8515 Miralani Drive
San Diego, CA 92126

argen.com | (800) 255-5524



Scan for Catalogs



Copyright © 2024 The Argen Corporation. All Rights Reserved.