

Ti Bases Gingival Heights

All platforms available in engaging and non-engaging

Series	Size	Gingival Heights
ASTRA TECH EV™* (SEV)	3.6	1.0mm, 3.0mm
	4.2	1.0mm, 3.0mm
	4.8	1.0mm, 3.0mm
ASTRA TECH OsseoSpeed TX™* (AO)	3.0 NP	1.0mm, 3.0mm
	3.5/4.0 RP	1.1mm, 3.0mm
	4.5/5.0 WP	0.9mm, 3.0mm
BioHorizons®* (BIO)	3.5 NP	0.8mm
	4.5 RP	0.3mm
	5.7 WP	0.3mm
Biomet 3i Certain®* (BC)	3.4 NP	0.3mm, 1.5mm
	4.1 RP	0.3mm, 1.5mm
	5.0 WP	0.3mm, 1.5mm
Biomet 3i OSSEOTITE®* (BO)	3.4 NP	0.5mm
	4.1 RP	0.5mm
	5.0 WP	0.5mm
CAMLOG®* (CAM)	3.3	0.5mm
	3.8	0.3mm, 1.5mm
	4.3	0.3mm, 1.5mm
	5.0	0.3mm, 1.5mm
Dentsply Friadent XiVE®* (DX)	3.4 NP	0.5mm
	3.8 RP	0.5mm
	4.5 WP	0.5mm
	5.5	0.5mm
Hiossen ET/HG®* (HIO)	3.5	1.0mm, 3.0mm
	4.0-6.0	1.0mm, 3.0mm
Megagen AnyRidge®* (MA)	Universal	2.0mm, 3.0mm
Neodent® Grand Morse™* (NGM)	Universal	0.8mm, 1.2mm
NobelActive™* (NC)	3.5 NP	1.2mm, 3.0mm
	4.3/5.0 RP	1.2mm, 3.0mm
	5.5 WP	1.5mm

Series	Size	Gingival Heights
Nobel Branemark®* (NB)	3.5 NP	0.5mm
	4.1 RP	0.5mm
	5.1 WP	0.5mm
NobelReplace Select®* (NR)	3.5 NP	0.3mm
	4.3 RP	0.3mm
	5.0 WP	0.3mm
Straumann®* Bone Level (SB)	3.3 NC	1.2mm, 3.0mm
	4.1/4.8 RC	1.2mm, 3.0mm
Straumann® SynOcta®* (SS)	3.5 NN	0.5mm
	4.8 RN	0.4mm
	6.5 WN	0.3mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5 NP	0.8mm, 1.5mm
	4.5 RP	0.3mm, 1.5mm
	5.7 WP	0.3mm, 1.5mm

*Not a trademark of The Argen Corporation.

The Argen Corporation

8515 Miralani Drive, San Diego CA 92126
argen.com | customersuccess@argen.com | (800) 255-5524



MKT-18953-Ver3-01/24