

Multi-Unit Abutment Gingival Heights and Angles

Series	Size	Multi-Unit Platform	Angle	Gingival Heights
Astra Tech OsseoSpeed® TX* (AO)	3.5/4.0	4.8	Straight	1.5 mm, 2.5 mm
	4.5/5.0	4.8	Straight	1.5 mm, 2.5 mm
Neodent® Grand Morse™*	Universal	4.8	Straight	1.5 mm, 2.5 mm
	Universal	4.8	17°	2.5 mm, 3.5 mm
	Universal	4.8	30°	3.5 mm, 4.5 mm
NobelActive®* (NC)	3.5 NP	4.8	Straight	1.5 mm, 2.5 mm
	3.5 NP	4.8	17°	2.5 mm, 3.5 mm
	3.5 NP	4.8	30°	3.5 mm, 4.5 mm
	4.3/5.0 RP	4.8	Straight	1.5 mm, 2.5 mm
	4.3/5.0 RP	4.8	17°	2.5 mm, 3.5 mm
	4.3/5.0 RP	4.8	30°	3.5 mm, 4.5 mm
NobelReplace® Select™* (NR)	3.5 NP	4.8	Straight	1.0 mm, 2.0 mm
	4.3 RP	4.8	Straight	1.0 mm, 2.0 mm
	5.0 WP	4.8	Straight	1.0 mm, 2.0 mm
Straumann® Bone Level* (SB)	3.3 NC	4.8	Straight	1.5 mm, 2.5 mm
	3.3 NC	4.8	17°	2.5 mm, 3.5 mm
	3.3 NC	4.8	30°	3.5 mm, 4.5 mm
	4.1/4.8 RC	4.8	Straight	1.5 mm, 2.5 mm
	4.1/4.8 RC	4.8	17°	2.5 mm, 3.5 mm
	4.1/4.8 RC	4.8	30°	3.5 mm, 4.5 mm
Zimmer Tapered Screw-Vent®* (ZS)	3.5	4.8	Straight	1.0 mm, 2.0 mm
	3.5	4.8	17°	2.5 mm, 3.5 mm
	3.5	4.8	30°	3.5 mm, 4.5 mm
	4.5	4.8	Straight	1.0 mm, 2.0 mm
	4.5	4.8	17°	2.5 mm, 3.5 mm
	4.5	4.8	30°	3.5 mm, 4.5 mm

*Not a trademark of The Argen Corporation.