

3D PRINTED DENTURES

A New Reason to Smile.

ARGEN[®] Digital Dentures

Argen 3D printed digital dentures and try-ins use the latest additive manufacturing technology and FDA-cleared materials. The CAD/CAM process creates true-to-impression quality restorations that improve accuracy for a better fit.

Argen offers easy file upload and quick turnaround times for accelerated production.



Assembled Digital Denture

Includes the printed pink base and printed tooth shade arch that has been assembled, surface treated, and UV cured. The denture will need to be polished and high shined at the lab.

Base shades: Original Pink, Dark Pink, Light Pink, Reddish Pink

Tooth shades: A1, A2, A3, A3.5, B1, B2



Try-In Digital Denture

Includes the printed tooth material in a 1-piece, monolithic denture. May require polishing at the lab.

Color: White



Place your order online
at **argen.com**

