

THE MOST ADVANCED ZIRCONIA TECHNOLOGY AVAILABLE

ARGEN  ANTERIOR
Super Translucent Zirconia



Fulfill all Prescriptions with 6 Value Shade Discs

VALUE SHADE DISC	DESIRED SHADE
Value Disc 1	A1 B1
Value Disc 2	A2 B2 A3
Value Disc 3	C1 D2
Value Disc 4	C2 D3 D4
Value Disc 5	A3.5 B3 B4
Value Disc 6	A4 C3 C4

- Exciting Alternative to Lithium Disilicate
- Up to 3 unit Anterior Bridges
- 50% Translucency
- 765 MPa Strength

MADE IN THE
USA

Place your order online at
argen.com

 **ARGEN**[®]
(800) 255-5524