



ARGEN^{DS} Denture Systems

ArgenDS HI-Base Denture Milling Discs

ArgenDS High Impact denture discs are used to mill premium dentures. Our HI-Base denture discs offer a highly esthetic, fracture resistant solution for your premium denture base needs. These strong, vibrant discs require minimal finishing and polishing.

ArgenDS HI-Base milling discs are 510(k) cleared and are compatible with 4-axis and 5-axis milling machines such as VHF, Roland, Imes-Icore, and others.

DIAMETER	THICKNESSES	SHADES
98 mm	25 mm 30 mm	Original, Light Reddish Pink

MATERIAL PROPERTIES	
ISO 20795 Ultimate Flexural Strength	90 MPa
ISO 20795 Flexural Modulus	2500 MPa
Residual monomer	<1% by weight
Water sorption	17 µg/mm ³
Water solubility	0.6 µg/mm ³



argen.com | (800) 255-5524



MKT-18477-Ver1-05/20