

Now Available

Argen-Milled IPS e.max

Outsource your milled IPS e.max cases to Argen.

Why mill with Argen?

Argen's state-of-the-art manufacturing facility and fast turn-around time allows you to benefit from an end-to-end digital workflow while bypassing cost and time barriers. Argen has the ability to produce the widest range of quality digital products that now includes Argen-Milled IPS e.max. Outsourcing your cases to Argen also gives you full access to our industry-leading technical support team. Relax and let Argen mill your e.max.

What's available?

- 4 Translucency Levels: HT, LT, MT, MO
- High Translucency Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, BL1, BL2, BL3, BL4, C1, C2, C3, C4, D2, D3, D4
- Low Translucency Shades: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, BL1, BL2, BL3, BL4, C1, C2, C3, C4, D2, D3, D4
- Medium Translucency Shades: A1, A2, A3, B1, BL2, BL3, BL4
- Medium Opacity Shades: 0, 1, 2, 3, 4
- **SINGLE UNITS ONLY**
- 2-day turn-around

What to expect?

Lab receives unit in blue state, uncrystallized.



As received
in-lab



After
crystallization



After a
quick glaze



(800) 255-5524

argen.com