

ARGEN[®]DS Denture Systems

ArgenDS HI-Base Denture Milling Discs

ArgenDS High Impact denture discs are used to mill premium dentures. Our HI-Base denture discs offer a highly esthetic, fracture resistant solution for your premium denture base needs. These strong, vibrant discs require minimal finishing and polishing.

ArgenDS HI-Base milling discs are 510(k) cleared and are compatible with 4-axis and 5-axis milling machines such as VHF, Roland, Imes-Icore, and others.

DIAMETER	THICKNESSES	SHADES
98 mm	25 30 mm	Original, Light Reddish Pink



ArgenDS Digital Denture

Argen 3D printed digital dentures and try-ins use the latest additive manufacturing technology and FDA-cleared materials. The CAD/CAM process creates true-to impression quality restorations that improve accuracy for a better fit.

Assembled Denture

Includes the printed pink base and printed tooth shade arch that has been assembled, surface treated, and UV cured. The denture will need to be polished and high-shined at the lab.

PINK BASE SHADES

Original Pink
Dark Pink
Light Pink
Reddish Pink

TOOTH SHADES

A1, A2, A3, A3.5, B1, B2



Try-In Denture

Includes the printed tooth material in a one-piece, monolithic denture. May require polishing at the lab.

COLOR

White



argen.com | (800) 255-5524



MKT-18434-Ver1-02/20