

# Add-Ons and Additional Seats

---

Add-Ons	Additional DSP or C&B Seat Needed
Full Denture 2 <sup>nd</sup> Seat	×
RPD 2 <sup>nd</sup> Seat	×
Model Builder Trios 2 <sup>nd</sup> Seat	×
Model Builder Allscans 2 <sup>nd</sup> Seat	×
Smile Designer 2 <sup>nd</sup> Seat	×
Splint Studio 2 <sup>nd</sup> Seat	×
Abutment Designer 2 <sup>nd</sup> Seat	✓
Implant Bar and Bridge Designer 2 <sup>nd</sup> Seat	✓

- Dentures and RPDs use a full denture library that is dependent on a denture tooth library from denture teeth companies. These libraries are different from the Crown and Bridge libraries.
- Model Builder and Splint Studio do not use a smile library therefore they can be open at the same time.
- Abutment and Implant Bar and Bridge use the dental designer smile library which requires an additional platform seat.