

ARGEN[®] IS

A better abutment.

Dentists and dental laboratories trust Argen to manufacture the highest quality custom prosthetic solutions compatible with the leading implant systems. Argen's state-of-the-art manufacturing facility in San Diego, CA provides end-to-end manufacturing allowing us to maintain complete control of the process. Our CNC Swiss Screw machines mill abutment blanks and screws from Grade 23 titanium, which is the highest purity medical grade titanium alloy available. All components are rigorously inspected using a Nikon Vision System that validates their specifications and tolerances. We guarantee the quality and fit of our abutments and offer an extensive warranty that covers both the implant and the abutment.

All ArgenIS blanks and screws are FDA 510(k) cleared. This clearance means Argen's implant components are certified compatible with and equivalent to the OEM (Original Equipment Manufacturer) parts. Dentists can rest assured that ArgenIS custom abutments meet the same quality, tolerance standards, and compatibility requirements of OEM units.



Available from your dental laboratory partner.



End-to-End
Manufacturing



FDA 510(k) Cleared /
ISO 13485



Advanced Milling
Strategies



Robotic
Automation



Extensive
Warranty



Watch Argen's ISA
manufacturing video.

Compatible with
16 of the leading
implant systems.



Why doctors choose Argen's Implant Specific Abutments.

	Argen's ISA* Abutments	Various OEM Abutments
FDA 510(k) Cleared	✓	✓
MDSAP Certified	✓	Varies
ISO 13485	✓	Varies
Warranty (Implant & Abutment)	✓	Varies
Strict Tolerances	✓	✓
True-to-Design Margins	✓	✓
Split Files	✓	Varies
FDA 510(k) Cleared 2 Final Screws	✓	Varies
Better Value	✓	—

*Implant Specific Abutment



Ask your dental laboratory for
Argen's Implant Specific Abutments



MKT-18080-Ver3-06/24