

Argen Titanium Intraoral Scan Bodies

Frequently Asked Questions

Q) What material are the scan bodies made out of?

- They are made of Ti 6AL 4V-ELI TITANIUM ALLOY (Grade 23)

Q) Do the scan bodies come with a screw?

- Yes, they come with a captured Ti screw

Q) Are the scan bodies coated?

- Yes, they have a gold colored ZrN coating to enhance scan ability and durability

Q) Are the scan bodies laser etched?

- Yes, they are marked with the series and size of the corresponding implant

Q) Can these scan bodies be used by the dentist?

- Yes, these scan bodies can be used by both the lab and the clinician

Q) Are the CAD files in the 3Shape and Exocad libraries?

- Yes, the CAD files are identical to the Gen 2 PEEK scan bodies, so they will work with the same libraries

Q) Do I need a new DME?

- No, Since the CAD files are the same as PEEK, the existing DME will work

Q) Are these scan bodies reusable?

- Yes, they are reusable and must be sterilized in between uses

Q) How many times can I sterilize the units?

- Consult the IFU for guidance

Q) Does Argen have a recommended sterilization cycle?

- Yes, consult the IFU

Q) Are there any contraindications to using Argen's Titanium Scan Bodies?

- Patients who are medically unfit for an oral surgical procedure
- Patients who are allergic or hypersensitive to titanium alloy Ti 6Al-4V ELI (Grade 23)

Q) What is the sales price?

- Contact your sales rep

Argen Titanium Intraoral Scan Bodies

Frequently Asked Questions

Q) Is there an IFU and SDS available?

- Yes, both are available on Argen's website, argen.com

Q) How can I get these new scan bodies?

- They are available on argen.com

Q) What Series and Platforms are Argen Titanium Scan Bodies available for?

- Astra Tech EV*
- Astra Tech OsseoSpeed TX®*
- Biomet 3i Certain®*
- BioHorizons®*
- CAMLOG®*
- Hiossen ET/HG*
- MegaGen AnyRidge®*
- NobelActive®*
- Neodent® Grand Morse™*
- NobelReplace Select®*
- Straumann® Bone Level*
- Straumann® synOcta® Tissue Level*
- Zimmer Tapered Screw-Vent®*

Q) Do Argen Titanium Scan Bodies come sterile?

- No, consult the IFU for guidance

*Not a trademark of The Argen Corporation.