

ArgenIS Implant System and Scan Bodies Compatibility

Implant Info			Indication Compatibility				Scan Body Compatibility						
System	Code	Size	Ti	Hybrid	ASC	Multi-Unit Abutment	ArgenIS **	Multi-Unit ***	NT **	3Shape	Elos	Elos IO	Dess IO
Astra Tech EV	SEV	3.0, 5.4	X	X			X		X	X		X	X
	SEV	3.6, 4.2, 4.8	X	X	X		X		X	X		X	X
Astra Tech OsseoSpeed® TX*	AO	3.0 NP		X									X
	AO	3.5/4.0 RP, 4.5/5.0 WP	X	X	X	X	X	X	X	X	X	X	X
BioHorizons®*	BIO	3.0	X				X		X	X			X
	BIO	3.5 NP, 4.5 RP, 5.7 WP	X	X			X		X	X	X	X	X
Biomet 3i Certain®**	BC	3.4 NP, 4.1 RP, 5.0 WP	X	X	X		X		X	X	X	X	X
Biomet 3i Osseotite®**	BO	3.4 NP, 5.0 WP		X	X				X	X	X	X	X
	BO	4.1 RP	X	X	X				X	X	X	X	X
CAMLOG®*	CAM	3.3, 5.0	X	X			X		X	X	X	X	X
	CAM	3.8, 4.3	X	X	X		X		X	X	X	X	X
Dentsply Friadent XiVE®**	DX	3.4 NP	X		X				X	X	X		X
	DX	3.8 RP, 4.5 WP	X		X				X	X	X	X	X
	DX	5.5	X						X	X	X	X	X
Hiossen ET/HG*	HIO	3.5, 4.0/4.5/5.0/6.0	X	X			X		X				X
MegaGen AnyRidge®**	MA	Universal	X	X	X		X						X
Neodent® Grand Morse™**	NGM	Universal	X	X	X	X	X	X					X
NobelActive®**	NC	3.5 NP, 4.3/5.0 RP	X	X	X	X	X	X	X	X	X	X	X
	NC	5.5 WP	X	X			X			X			X
Nobel Brånemark System®**	NB	3.5 NP, 4.1 RP	X	X	X				X	X	X	X	X
	NB	5.1 WP	X	X					X	X	X	X	X
NobelReplace®**	NR	3.5 NP, 4.3 RP, 5.0 WP	X	X	X	X	X	X	X	X	X	X	X
	NR	6.0		X			X		X	X			X
Straumann® Bone Level*	SB	3.3 NC, 4.1/4.8 RC	X	X	X	X	X	X	X	X	X	X	X
Straumann® synOcta®**	SS	3.5 NN		X					X	X			X
	SS	4.8 RN, 6.5 WN	X	X	X		X		X	X	X	X	X
Zimmer Tapered Screw-Vent®*	ZS	3.5 NP, 4.5 RP	X	X	X	X	X	X	X	X	X	X	X
	ZS	5.7 WP	X	X	X		X		X	X	X	X	X

* Not trademarks of The Argen Corporation. ** Sold by Argen. *** Sold by Argen, one size fits all.