

# ArgenIS Angled Screw Channel Abutments

## Screw Torque Setting

The following suggested screw torque settings are for reference and have been provided for your convenience.

Implant System	Series	Torque Setting
ASTRA TECH EV™*	SEV Series	25 N-cm
	SEV Series 4.8	25 N-cm
ASTRA TECH OsseoSpeed TX™	AO Series 3.4	15 N-cm
	AO Series 3.8	20 N-cm
	AO Series 4.5	25 N-cm
Biomet 3i Certain®*	BC Series	20 N-cm
Biomet 3i OSSEOTITE®*	BO series	20 N-cm
CAMLOG®*	CAM Series	20 N-cm
Dentsply Friadent XiVE®*	DX Series	24 N-cm
Megagen AnyRidge®*	MA Series	35 N-cm
NobelActive®*	NC Series	35 N-cm
NobelReplace Select®*	NR Series	35 N-cm
Nobel Branemark®*	NB Series	35 N-cm
Straumann Bone Level®*	SB Series	35 N-cm
Straumann SynOcta®*	SS Series	35 N-cm
Zimmer Tapered Screw-Vent®*	ZS Series	30 N-cm

\* Not trademarks of the Argen Corporation.

### ARGEN EXTENSIVE ABUTMENT WARRANTY\*

Argen Digital dental abutments are backed by a Limited Warranty for titanium custom abutments and zirconia hybrid abutments. In addition, if an implant supplier does not honor their Warranty due to the use of an Argen abutment, Argen will replace both the abutment and the implant.

*\*See complete warranty details at [www.argen.com](http://www.argen.com)*