

ArgenIS Screw Torque Settings

Implant Systems	ArgenIS Screws					
	Series	Titanium Custom Abutment	Hybrid Custom Abutment	Angled Hybrid Abutment	Multi-Unit Abutment	Multi-Unit Abutment Ti Base Cap
ASTRA TECH EV™*	SEV	25 Ncm	25 Ncm	25 Ncm	N/A	N/A
ASTRA TECH OsseoSpeed TX™*	AO	25 Ncm	15 Ncm (3.0) 20 Ncm (3.5/4.0) 25 Ncm (4.5/5.0)	20 Ncm (3.5/4.0) 25 Ncm (4.5/5.0)	15 Ncm	15 Ncm
BioHorizons*	BIO	30 Ncm	30 Ncm	N/A	N/A	N/A
Biomet 3i Certain®*	BC	20 Ncm	20 Ncm	20 Ncm	N/A	N/A
Biomet 3i OSSEOTITE®*	BO	35 Ncm	20 Ncm	20 Ncm	N/A	N/A
CAMLOG®*	CAM	20 Ncm	20 Ncm	20 Ncm	N/A	N/A
Dentsply Friadent XiVE®*	DX	25 Ncm	24 Ncm	24 Ncm	N/A	N/A
Hiossen ET / HG*	HIO	20 Ncm (3.5) 30 Ncm (4.0-6.0)	20 Ncm (3.5) 30 Ncm (4.0-6.0)	N/A	N/A	N/A
Megagen AnyRidge*	MA	N/A	35 Ncm	N/A	N/A	N/A
NobelActive®*	NC	20 Ncm (3.5) 35 Ncm (4.3-5.5)	35 Ncm (3.5-5.5)	35 Ncm	35 Ncm	15 Ncm
NobelReplace Select®*	NR	35 Ncm	35 Ncm	35 Ncm	35 Ncm	15 Ncm
Nobel Branemark®*	NB	35 Ncm	35 Ncm	35 Ncm	N/A	N/A
Straumann Bone Level®*	SB	35 Ncm	35 Ncm	35 Ncm	35 Ncm	15 Ncm
Straumann SynOcta®*	SS	35 Ncm	35 Ncm	35 Ncm	35 Ncm	15 Ncm
Zimmer Tapered Screw-Vent®*	ZS	30 Ncm	30 Ncm	30 Ncm	30 Ncm	15 Ncm

*Not trademarks of the Argen Corporation.