

3Shape Scanner PC Options



Standard PC

Recommended for E1/E2

Lenovo

- I5 - 9400
- Nvidia Quadro P620 2GB
- 512GB M2 SSD
- 16GB RAM



High Performance PC

Recommended for E3/E4/D2000

Lenovo

- I7 - 9700
- Nvidia Quadro P1000 4GB
- 512GB M2 SSD
- 32GB RAM



Ultra Performance PC

Lenovo

- I7 - 9700
- Nvidia Quadro P2000 5GB
- 1024GB M2 SSD
- 64GB RAM