


STRONG. TRANSLUCENT. LAYERED.

More sizes, accurate shades, incredible strength.



ARGEN Z HT+
HT+ MULTILAYER

ArgenZ HT+ Multilayer pairs our proprietary layer blending technology with the strength and consistent shade match of our proven high translucent plus material. The seamless blend of this homogenous 4Y zirconia allows you to create large cases that mirror the natural gradient of dentin and enamel without visible layer lines or nesting challenges. Suitable choice for single units to full arch restorations.

Glaze & Go!

Available in light shades
and high definition milling

TRANSLUCENCY
45%

STRENGTH
1250 MPa

19 SHADES
16 VITA Classical
+ 3 Bleach, 5 Light Shades*

7 THICKNESSES
12, 14, 16, 18,
20, 25, 30mm

DIAMETER
98mm

*The "L" or "Light" shades were formulated for implant-retained full-arch restorations.



argen.com



MKT-18551-Ver4-09/22