

# What's Hot with ArgenIS Implant Systems

**ARGEN<sup>IS</sup>**  
Implant Systems



## New Ti Bases

- Engaging and non-engaging
- Multiple gingival heights
- Straight and angled screw channel
- Textured surface for strong adhesion to the restoration
- New platform: Neodent<sup>®</sup> Grand Morse<sup>™</sup>

## New Lab Analogs

- Color anodized for ease of size identification and to match OEM implant parts

## Temporary abutments

- Horizontal grooves allow for interlocking with multiple restoration materials
- Engaging and non-engaging

## Transfer Coping

- Can be used for open or closed tray impressions

## New Multi-Unit System

- Multi-Unit Abutments
  - Straight and angled
  - Multiple gingival heights
- Multi-unit Ti Bases
  - Straight offered in 6 mm and 8 mm chimney heights
  - Angled offered in one chimney height
  - RP size
- Multi-Unit Temporary Abutments
- Multi-Unit Transfer Coping
  - Can be used for open or closed tray impressions
- Multi-Unit analogs
- Multi-Unit scan bodies
- Multi-Unit key, adapter for multi-unit key, and driver

## NobelActive Titanium Abutment - Size 5.5

- NobelActive titanium abutment now available in size 5.5

## ArgenIS Angled Screw Channel

- Access hole can be placed anywhere between 0° to 25° in a 360° radius
- Engaging and non-engaging
- Features an engineered surface coating for optimal cement retention



[argen.com](http://argen.com) | (800) 255-5524



MKT-188876-Ver3-01/24